BYPASS FUME HOODS

HMF-BG

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Real-Time Face Velocity Monitoring

Continuously tracks and displays face velocity in real-time using high-precision sensors. If velocity drops below the safe threshold, the system immediately alerts the user to ensure proper sash operation and maintain lab safety.

Example: "0.00 m/s", "Low face velocity. Close the sash"

Intuitive Interface with High Readability

Key parameters such as velocity, time, and efficiency are presented in large digits and visual icons on a high-resolution TFT LCD. This enhances operational clarity and reduces the risk of user error.

Integrated Equipment Control

Operators can control key fume hood functions directly via the touchscreen:

- Lamp ON/OFF
- Blower ON/OFF
- Door open/close status monitoring
- Emergency response activation

Comprehensive Environmental Data Monitoring

In addition to airflow, the system monitors various environmental metrics:

- Temperature & humidity levels
- VOCs (Volatile Organic Compounds) concentration
- Cumulative operation time and system efficiency

Safety-First Visual Alert System

Dynamic alert banners with color-coded warnings notify users of hazardous conditions. Immediate visual feedback minimizes response time and reinforces safety protocols.

Mobile-Controlled Operation with Bluetooth Connectivity

The Huco touchscreen system is seamlessly integrated with a dedicated mobile application, enabling real-time monitoring and control of core functions via smartphone. Its Bluetooth-enabled wireless communication allows for convenient remote access without complex network configuration, enhancing operational flexibility and user convenience in dynamic industrial settings.