



# Viscometry

VISCOSYSTEM® | AVS® | VISCOCLOCK | VISCOMETERS

SI Analytics

a **xylem** brand

# Capillary viscometry from Xylem - know-how from the very beginning



## Innovative capillary viscometry - from the outset

The viscosity of Newtonian fluids can be most precisely determined using capillary viscometers. This method of measurement, measures the time taken for a defined quantity of fluid to flow through a capillary with a known diameter and known length. With the industrial production of such precisely calibrated capillary viscometers, we have created the conditions to enable this measuring method to establish itself worldwide as a reliable procedure.

With the development of the first automatic measuring systems, we replaced the stopwatch with automatic registration of the fluid at the start of the 1970's.

To rule out systematic errors in automatic meniscus detection, viscometers are available which are specially calibrated for this type of automatic cycle time measurement and therefore comply with the relevant standard specifications.

Our capillary viscometers are the worldwide basis for precise viscosity measurements of Newtonian fluids.

In addition to automatic time measurement, AVS® measuring systems have further automations that make capillary viscometry much easier:

- Pneumatic pumping of the liquid in the viscometer, so that automatically repeated measurements can be carried out, from which an average value and the final result are calculated
- Waste system: Automatic emptying and flushing of viscometers
- The automatic filling of the viscometers on the AVS®Pro III automatic sampler

In addition to the accuracy of the measurement, all automated devices focus on user safety and the robustness of the system.

Further developments and improvements of viscometers, measuring instruments and accessories led to a range of products whose excellent performance is universally recognized. It is therefore no wonder that our viscosity measurement systems have become indispensable production control and quality insurance tools worldwide, whether in the mineral oil industry, for polymer manufacturers and processors, in the pharmaceutical or food industry.

# Content

|   | Page |
|---|------|
| <b>1. Measurement Devices</b>   |      |
| 1.1 Applications of the AVS® measurement systems                              | 4    |
| 1.2 ViscoClock plus   | 6    |
| 1.3 ViscoPump III - The core of all AVS® systems                              | 10   |
| 1.4 AVS® 470  | 12   |
| 1.5 AVS® 370  | 16   |
| 1.6 Software WinVisco 4   | 22   |
| 1.7 AVS® Pro III  | 26   |
| 1.8 CT 72 Thermostat Series - Transparent Thermostats conforming to Standards | 34   |
| <b>2. Viscometer</b>  |      |
| 2.1 Viscometers and their range of use  | 38   |
| 2.2 Ubbelohde viscometers, normal form (DIN)                                  | 39   |
| 2.3 Ubbelohde viscometers, normal form (ASTM)                                 | 40   |
| 2.4 Ubbelohde viscometers,<br>with additional tube and threads                | 41   |
| 2.5 Ubbelohde viscometers with TC sensors                                     | 42   |
| 2.6 Micro-Ubbelohde viscometers with TC sensors                               | 44   |
| 2.7 Micro-Ubbelohde viscometers (DIN)   | 45   |
| 2.8 Cannon-Fenske viscometers   | 46   |
| 2.9 Ostwald Viscometer  | 48   |
| <b>3. Accessories</b>   |      |
| 3.1 Accessories for viscometers   | 49   |
| 3.1.1 Brackets and stands   | 49   |
| 3.1.2 Temperature stabilization jackets                                       | 51   |
| 3.1.3 LabPump   | 52   |
| 3.2 Accessories for measurement devices                                       | 52   |
| 3.2.1 Control thermometers  | 52   |
| 3.2.2 Safety devices  | 54   |
| 3.2.3 Filtration system ProClean II   | 56   |
| 3.2.4 Thermostat vessel for ViscoClock <i>plus</i>                            | 57   |
| 3.2.5 Burette Titronic® 500 for sample preparation                            | 58   |
| 3.2.6 AVS® measuring stands   | 60   |

# 1 Measurement Devices

## 1.1 Applications AVS® systems

### AVS® measurement systems within quality assurance systems

| Business sector                        | Product   | Example  |
|--|---|--|
| Automotive engineering                 | motor oil (fresh and used)<br>uncured plastics                                      | light weight construction parts                            |
| Brewery                                | original wort<br>hop-wort   | beer<br>beer   |
| Electrical engineering and electronics | uncured plastics of all types   | chips, casings   |
| Power supply                           | turbine oil<br>transformer oil  | generators   |
| Plastics manufacturers                 | uncured plastics of all types   | Polyamid (Nylon)   |
| Plastics processors                    | uncured plastics of all types   | injection moulding of polyester and polyamide              |
| Food industry                          | starch  | instant flour thickeners                                   |
|  | gelatin   | jelly bears  |
|  | packaging materials   | yogurt containers  |
|  | milk products   | yogurt drink   |
|  | fruit and fruit juice concentrates<br>gelatinizing agents                           | pectin   |
| Aviation                               | uncured plastics of all types<br>fuels<br>hydraulic fluids                          | kerosene<br>horizontal stabilizers and undercarriages      |
| Mechanical engineering                 | mold oil<br>hardening emulsions<br>hydraulic fluids                                 | mill trains<br>stamp shops                                 |
| Medicine                               | body fluids<br>hyaluronic acid<br>tinctures and drops<br>blood substitute materials | blood, bile<br>Na-Hyaluronat<br>nose, eyes<br>blood plasma |
| Mineral oil                            | light motor oil<br>turbine oil<br>liquid fuels of all types                         | gasoline, diesel fuel, kerosene (jet fuel)                 |
| Textile                                | uncured plastics of all types<br>cotton   | for mixed fibers   |
| Entertainment                          | uncured plastics  | CDs, DVDs made of polycarbonate                            |

The table on the right illustrates the extensive area of high polymer plastics and the large variety of testing methods.

# Polymer applications for the AVS® measurement systems

## Recommended capillary sizes and AVS® measurement systems

For polymer analytics, mostly Ubbelohde viscometers are applied (type with suspended level). In general, viscometers according to DIN 53 000 as well as ASTM D446 can be used, as there is only a slight difference between both designs. Especially ISO standards for polymer analytic generally allow both designs. However, in case of ASTM standards, e.g. ASTM D4603, the ASTM viscometer design is more common.

| Type                            | Abbr. | Solvent                                | Capillary size         |         | Temperature | Operating Standards        | Suitability of the AVS® measurement systems |     |     |         |
|---------------------------------|-------|--|------------------------|---------|-------------|----------------------------|---|-----|-----|---------|
|                                 |       |  | DIN                    | ASTM    |             |                            | VC*   | 370 | 470 | Pro III |
| Cellulose                       | CI    | EWN                                    |                        |         | 20 °C       | SNV 195598                 | ■   | ■   | ■   | ■       |
|                                 |       | Cuen (CED)                             |                        |         | 20 °C       | DIN EN 60450               |   |     |     |         |
|                                 |       | Cuen (CED)                             |                        |         | 20 °C       | ASTM D 4243                |   |     |     |         |
|                                 |       | Cuen (CED)                             |                        |         | 25 °C       | ISO 5351                   |   |     |     |         |
|                                 |       | Cuen (CED)                             |                        |         | 25 °C       | ASTM; D 1795               |   |     |     |         |
|                                 |       | Cuen (CED)                             |                        |         | 25 °C       | SCAN CM 15:99              |   |     |     |         |
|                                 |       | Cuen (CED)                             |                        |         | 25 °C       | TAPPI T230-0M99            |   |     |     |         |
| Cellulose acetate               | CA    | Dichlormethane/methanol (90:10)        | 0c<br>I<br>I Micro     | 0C<br>1 | 25 °C       | ASTM D817                  | ■   | ■   | ■   | ■       |
| Polyamide                       | PA    | Sulphuric acid (96%)                   | II<br>IIc              | 2<br>2C | 25 °C       | ISO 307                    | ■   | ■   | ■   | ■       |
| Polyamide                       | PA    | Formic acid (90%)                      | I<br>Ic                | 1<br>1C | 25 °C       | ISO 307                    | ■   | ■   | ■   | ■       |
| Polyamide                       | PA    | m-cresol                               | II<br>IIc              | 2<br>2C | 25 °C       | ISO 307                    | ■   | ■   | ■   | ■       |
| Polybutylene terephthalate      | PBT   | Phenol/dichloro benzene (50:50)        | Ic<br>II               | 1C<br>2 | 25 °C       | ISO 1628-5                 | ■   | ■   | ■   | ■       |
| Polycarbonate                   | PC    | Dichloromethane                        | 0c<br>I                | 0C<br>1 | 25 °C       | ISO 1628-4                 | ■   | ■   | ■   | ■       |
| Polyethylene                    | PE    | Decahydro-naphthalene                  | 0a<br>I                | 0B<br>1 | 135 °C      | ISO 1628-3<br>ASTM D 1601  | ■   | ■   | ■   | ■       |
| Polyethylene terephthalate      | PET   | Phenol/1,1,2,2-Tetrachlorethan (60:40) | Ia                     | 1B      | 25 °C       | ASTM D 4603                | ■   | ■   | ■   | ■       |
| Polyethylene terephthalate      | PET   | m-cresol                               | II<br>IIc<br>IIc Micro | 2<br>2C | 25 °C       | ISO 1628-5                 | ■   | ■   | ■   | ■       |
| Polyethylene terephthalate      | PET   | Phenol/dichloro benzene (50:50)        | Ic<br>II               | 1C<br>2 | 25 °C       | ISO 1628-5                 | ■   | ■   | ■   | ■       |
| Polyethylene terephthalate      | PET   | Dichloroacetic acid                    | Ia<br>II<br>IIc Micro  | 1B      | 25 °C       | ISO 1628-5                 | ■   | ■   | ■   | ■       |
| Polymethyl methacrylate         | PMMA  | Chloroform                             | 0c<br>I Micro          | 0C      | 25 °C       | ISO 1628-6                 | ■   | ■   | ■   | ■       |
| Polypropylene                   | PP    | Decahydro-naphthalene                  | 0a<br>I                | 0B<br>1 | 135 °C      | ISO 1628-3                 | ■   | ■   | ■   | ■       |
| Polystyrene                     | PS    | Toluene                                | I<br>Ic                | 1<br>1c | 25 °C       |                            | ■   | ■   | ■   | ■       |
| Polysulphone                    | PSU   | Chloroform                             | 0c                     | 0C      | 25 °C       |                            | ■   | ■   | ■   | ■       |
| Polyvinyl chloride              | PVC   | Cyclohexanone                          | Ic                     | 1C      | 25 °C       | ISO 1628-2,<br>ASTM D 1243 | ■   | ■   | ■   | ■       |
| Styrene-acrylonitrile copolymer | SAN   | Ethyl methyl ketone                    | 0c<br>I                |         | 25 °C       |                            | ■   | ■   | ■   | ■       |
| Styrene-butadiene copolymer     | SB    | Toluene                                | 0c<br>I                |         | 25 °C       |                            | ■   | ■   | ■   | ■       |

■ excellent suitability   ■ can be used   ■ limited suitability for application related reasons

This table makes no claim to completeness.

VC\* = ViscoClock plus  
The highlighted capillary size is specified in the standard

# 1.2 ViscoClock *plus*

## Measurement *plus* data storage

The ViscoClock *plus* is an electronic timing unit for glass capillary viscometers used to determine kinematic and relative viscosity. Succeeding the well-proven ViscoClock, the new instrument features data storage and simpler handling. The ViscoClock *plus* is especially designed for Ubbelohde type viscometers which are well-known for highest precision.

### The ViscoClock *plus*

The ViscoClock *plus* automatically measures the flow time of temperature-stabilized liquids in capillary viscometers by means of infrared light barriers: the manual measurement with a stopwatch becomes obsolete.

The viscometer including a sample is inserted into the ViscoClock *plus* and immersed into a thermostatic bath for temperature stabilization. After thermostating, the sample is pumped into the measuring bulb, and the flow time is detected automatically. The large display enables easy read-off of flow times and additional information: date, time, sample ID and viscometer ID.



*ViscoClock plus - Head*



*The viscometer stand:  
high-performance  
plastic PPA*



*Well arranged:  
The new display*

## Automatic measurement of flow times

The ViscoClock *plus* is designed for SI Analytics® Ubbelohde, Micro Ubbelohde and Micro Ostwald viscometers. The flow time is measured automatically by two infrared light barriers which detect the passing liquid meniscus. The repeatability of the automatic time measurement is considerably higher in comparison to the measurement using a stop watch. Therefore some viscometry standards allow a flow time reduction in case of automatic flow time measurement.

## Properties and materials

The ViscoClock *plus* can be used for measuring temperatures ranging from -40 °C to 150 °C. The stand of the ViscoClock *plus* is made of high quality polymer PPA. For temperature stabilization in a thermostatic bath, the following liquids are suitable: Water, alcohol, glycol, paraffin oil, and silicon oil.

The electronic measuring unit is built-in to a PP casing.

## Easy handling

The ventilation of Ubbelohde viscometers is managed by an electromechanic valve which makes handling easier in comparison to the mechanical mechanism of the previous ViscoClock.

## Data storage

The measuring results of the ViscoClock *plus* can be stored on a USB flash drive including date, time and sample/viscometer ID. The data are stored as pdf (non-editable) and csv (editable). Alternatively, for data transfer the ViscoClock *plus* can be connected to a printer (TZ 3863) or a PC.



- ▶ Automatic and precise flow time measurement for a low price
- ▶ Suitable for SI Analytics® Ubbelohde, Micro Ubbelohde and Micro Ostwald viscometers
- ▶ Data storage incl. time, date, viscometer and sample ID
- ▶ Stand is made of high performance engineering plastic PPA and enables measuring temperatures up to 150°C
- ▶ Electromagnetic venting valve for convenient handling of Ubbelohde viscometer
- ▶ Compatible with all SI Analytics® thermostatic bath types

Benefits  
ViscoClock *plus*

# ViscoClock *plus* - The *plus* for your measurements

## Sample and viscometer identification

To allocate the stored measuring results, the user can enter 2-digit numbers to the ViscoClock *plus* before measurement. These IDs - together with date and time - ensure an unambiguous assignment of the flow times.

## Absolute viscosity

To determine absolute kinematic viscosities, calibrated viscometers have to be used. To guarantee best accuracy, viscometers which were calibrated by automatic measurement should be used. The constant of automatic calibration can be slightly different in comparison to manual calibration, as the level of the light barriers may not be identical to the position of timing marks.

## Relative viscosity

In the analytics of plastics, for evaluation the relative viscosity is calculated, and depending on this also viscosity number (VN), intrinsic viscosity (IV) or the K value according to Fikentscher. For determination of relative viscosities, calibrated as well as non-calibrated viscometers can be used. For evaluation, the calibration constant is not required in this case.

## Ordering Information

| Type No.                        | Order No. | Description   | Page  |
|---------------------------------|-----------|---|-------|
| ViscoClock <i>plus</i>          | 285417900 | Timing unit for capillary viscometer. Including power supply 100-230V and hand pump                           | 9     |
| ViscoClock <i>plus</i> M1, 230V | 285417910 | ViscoClock <i>plus</i> and acrylic glass thermostatic bath CT72/P (230V) for temperatures +10 °C ... +60 °C   | 9, 39 |
| ViscoClock <i>plus</i> M1, 115V | 285417920 | ViscoClock <i>plus</i> and acrylic glass thermostatic bath CT72/P (115V) for temperatures +10 °C ... +60 °C   | 9, 39 |
| ViscoClock <i>plus</i> M2, 230V | 285417930 | ViscoClock <i>plus</i> and glass panelled thermostatic bath CT72/2 (230V) for temperatures -40 °C ... +150 °C | 9, 39 |
| ViscoClock <i>plus</i> M2, 115V | 285417940 | ViscoClock <i>plus</i> and glass panelled thermostatic bath CT72/2 (115V) for temperatures -40 °C ... +150 °C | 9, 39 |
| Thermostat vessel               | 285424400 | Thermostat vessel ViscoClock <i>plus</i>  | 57    |

# Technical Data - ViscoClock *plus*

|                                       |  |
|---------------------------------------|--|
| <b>Measuring range - Time</b>         | up to 999.99 s; resolution 0.01 s  |
| <b>Accuracy of time measurement</b>   | $\pm 0.01$ s/ $\pm 1$ digit; however no more precise than 0.1 %;<br>indicated as measuring uncertainty with a confidence level of 95 %   |
| <b>Measuring range - viscosity</b>    | 0.35 to 10,000 mm <sup>2</sup> /s (cSt)<br>the absolute, kinematic viscosity is additionally dependent on the uncertainty of the numerical value of the viscometer constant and on the measuring conditions, in particular the measuring temperature.  |
| <b>Display</b>                        | LCD graphic display (FSTN) 128 x 64 pixel, 51x31 mm (w x h)<br>seconds indication with 2 decimal digits after the decimal point, resolution 0.01 s   |
| <b>Voltage supply</b>                 | DC + 9 V   |
| <b>Power supply</b>                   | in accordance to class of protection III<br>degree of protection for dust and humidity IP 50 in accordance with DIN 40 050<br>Universal power supply TZ 1858: 100-240 V, 50-60 Hz (9 V, 550 mA)<br>not suitable for use in areas subject to explosion hazards  |
| <b>Interfaces</b>                     | USB Host to connect USB flash drive or printer (TZ 3863)<br>USB OTG to connect (PC), printer (TZ 3863) or USB flash drive  |
| <b>Plug Connections</b>               | socket for low voltage connection: coaxial power connector, inner diameter 2.1 mm, plus pole at inner contact for connection of Universal power supply TZ 1858<br>Type A USB connector<br>Type B mini USB connector  |
| <b>Ambient Conditions</b>             | Ambient temperature +10 to +40 °C for storage and transport<br>Operating temperature stand: -40 to +150 °C<br>electronic measuring unit: +10 to +40 °C<br>Humidity in accordance with EN 61 010, Part 1<br>max. relative humidity 80% for temperatures up to 31 °C,<br>decreasing linearly to 50% of relative humidity at a temperature of 40 °C |
| <b>Housing</b>                        | Materials stand: polyphthalamide (PPA)<br>casing: polypropylene (PP)<br>gaskets: silicone<br>Dimensions ~515 x 90 x 30 mm (H x W x D)<br>Weight ~450 g (without viscometer)<br>power supply unit: ~220 g   |
| <b>Country of origin</b>              | Federal Republic of Germany  |
| <b>CE symbol</b>                      | In accordance with low voltage guideline 2014/35/EU<br>Test regulation EN 61 010-1:2011-07 for laboratory instruments in accordance with EMC regulation 2014/30/EU<br>Test regulation EN 61 326 Part1:2012<br>In accordance with RoHS regulation 2011/65/EU<br>Test regulation EN 50 581:2013-02<br>FCC Symbol                                   |
| <b>Viscometer types</b>               | Ubbelohde (DIN; ISO; ASTM; Micro), Micro-Ostwald, type SI Analytics®   |
| <b>Transparent thermostatic baths</b> | The ViscoClock <i>plus</i> can be used in all SI Analytics® bath types   |

## 1.3 ViscoPump III

The ViscoPump exists currently in the 3rd generation. It is included in all AVS® devices and has several central functions:

- The signals of automatic meniscus detection - from measuring stands in the case of optical detection or TC viscometers - are evaluated to determine the flow time.
- With a built-in pump, the liquid is pneumatically pumped up into the measuring bulb of the viscometer.
- The available working modes are either the classic pressure or the suction mode. The operation (by suction of the liquid in the capillary tube) has advantages in some applications, in particular in automated measuring stations with a waste system.
- The built-in micro-processor controls the pumping process via several adjustable parameters.
- In the case of a connected waste system, the discharging by an external pump is also controlled by the ViscoPump III.



## Ordering information

| Type No. | Order No. | Description                         |
|----------|-----------|-------------------------------------|
| VZ 8561  | 285424060 | ViscoPump III for optical detection |
| VZ 8562  | 285424070 | ViscoPump III for TC detection      |

# The core of all AVS® systems

The measurement control by the ViscoPump is part of the modular device concept of AVS® systems:

- In multi-channel systems, one ViscoPump III controls the measurement at one measuring position.
- When servicing, a ViscoPump III can be easily replaced by the user.

Despite a new microcontroller, the new ViscoPump III is functionally compatible with its predecessors ViscoPump and ViscoPump II: Existing ViscoPump modules from older generations can be replaced by the new ViscoPump III. Also a mixed assembly in multi-channel systems, e.g. with ViscoPump II and ViscoPump III, is possible.

The ViscoPump III is available in two versions, for optical or TC detection.

The ViscoPump III is usually delivered in a package together with other components as a complete AVS® system, and only ordered separately in case of replacement or upgrade to additional measuring positions.

- ▶ Optimized command sequences of ViscoPump III compared to ViscoPump II
- ▶ Compatible to earlier versions
- ▶ Update possible via internal USB port
- ▶ Modular concept
- ▶ In case of service, the ViscoPump III can be easily replaced by the user.
- ▶ Pumping mode reversible, between suction and pressure mode
- ▶ Proven mechanical components (pump, valve) of high durability

Benefits  
ViscoPump III



## 1.4 AVS® 470

### Precise Capillary Viscometry – Easy, Flexible and Independent of a PC

- ▲ Perfectly equipped for automatic viscosity measurements

The AVS® 470 is a measuring system that includes everything you need to take precise and reproducible measurements. All types of viscosity calculations with polymer solutions are integrated into the device. A small keyboard allows you to enter additional data. A serial printer can be used to conveniently document your results.

So, in a minimum of space, you can set up a measuring station equal in every way to complex measuring installations in terms of precision and reproducibility.


- ▲ Simple and updateable Modular Concept

The AVS® 470 is equipped with a ViscoPump III and therefore limited to one measuring station.

Due to its modular design, the AVS® 470 can be used with a ViscoPump III, either for optical or TC detection.

You can use your existing accessories such as thermostats, stands or flow-through coolers. Also, virtually all SI Analytics® capillary viscometers can be used.





Preferred "Suction" or "Pressure"?

Applications in comparison

|  | "Pressure" | "Suction" |
|--|------------|-----------|
| highly viscous samples e.g. oils, polymers | ■          | ■         |
| Solvents: highly volatile                  | ■          | -         |
| (examples) Dichloromethane                 | ■          | -         |
| Chloroform                                 | ■          | -         |
| Sulfuric acid                              | -          | ■         |
| Dichloroacetic acid                        | -          | ■         |
| Toluene                                    | ■          | ■         |
| Hexafluoro-isopropanol                     | ■          | ■         |
| m-cresol                                   | -          | ■         |
| Formic acid                                | -          | ■         |
| Phenol / dichlorobenzene                   | -          | ■         |
| Phenol / Tetrachloroethane                 | -          | ■         |
| Without waste system                       | ■          | -         |
| With waste system                          | -          | ■         |

- ▶ Automatic and highly precise measurements - independent of a PC
- ▶ "Suction" and "pressure" measurements with the same system
- ▶ Simple data input and parameterization via included mini-keyboard
- ▶ GLP documentation compliant when connected to an optional printer

Advantages  
AVS<sup>®</sup> 470

# AVS® 470 – Precise and Reliable

## Working with the AVS® 470 is easy

The desired measurement method can be pre-selected and started on the device. The entire measurement is automatic to eliminate subjective measurement errors. Once the set pre-heating time is reached, the desired number of measurements are taken. The status of the measurements is continuously displayed.

If required, individual parameters may be input via an included keyboard. A serial printer can be used to print measurement logs.

The connections are on the front panel of the device for easy control. Over-pumping and oversuction are prevented by the use of an optional capacitive sensor.

## Clear user guidance, clear status – even without PC!

The print-out shows everything you need for reliable documentation of your test.

```

No. 1 = 77.20s
No. 2 = 77.21s
No. 3 = 77.20s

=====

*****
*                               *
*   ViscoSystem AUS470         *
*       protocol               *
*                               *
*****

method : absolute

Id : 11
lot: SIM Test sample
usr: A. Eich

measurements [s]
No. 1 = 77.20*
No. 2 = 77.21*
No. 3 = 77.20*

delta%choice = 0.01%
pre temp. time = 0min

average      = 77.203s
stand. dev. = 0.006

constant = 0.029999996

AbsVisc=2.3161mm^2/s

temperature: 25.00 C
date:      05/12/2017
time:      09h 47m 27s
=====
    
```

Individually determined readings

Indication of method set

Designation of specimen

Readings used for evaluation

Set equalization time

Corrected average running time

Calculated Viscosity

Charge Number


User

Set maximum permissible deviation from average


Average of running times

Viscosimeter constant


Operating temperature, date and time at time of test




After switching on the AVS® 470 a self test is run and then an entry prompt appears.



The parameters can be set in the test mode. In case of a required blank value, it is saved in the system.



All setup parameters can be preset conveniently, e.g. pressure/suction, pumping speed, waiting time between two tests, etc.



The readings can be read off conveniently on the display regardless of whether or not a printer is connected.

# Technical data

|                                 |   |   |  |
|---------------------------------|---|---|--|
| Measuring range (time)          | 5 s up to 9,999.99 s; resolution 0.01 s                                 |   |  |
| Measuring range (viscosity)     | pressure:   | 0.35 to 1,800 mm <sup>2</sup> /s (cSt)                  |  |
|                                 | suction:  | 0.35 to ~5,000 mm <sup>2</sup> /s (cSt)                 |  |
| Measured parameter              | flow-through time [s]   |   |  |
| Time measuring accuracy         | ± 0.01 %  |   |  |
| Measured value display          | LC-Display  |   |  |
| Display accuracy                | ± 0.01 s, ± 1 Digit, but not exceeding 0.1%                             |   |  |
| Pumping pressure                | fully automatically controlled  |   |  |
|                                 | suction up to ~-160 mbar, pressure up to ~+160 mbar                     |   |  |
| Preselectable tempering period  | 0 to 20 min   |   |  |
| Connections                     | Pneumatic connections   | threaded connections for viscometers                    |  |
|                                 | Electrical connections  | circular connector with bayonet lock for viscometer     |  |
|                                 |   | 4-pin DIN socket for TC viscometer                      |  |
|                                 |   | 4-pin circular connector for capacitive sensor          |  |
|                                 |   | 7-pin circular connector for AVS® 26, with bayonet lock |  |
|                                 | RS232-C interface   | 9-pin for serial printer                                |  |
|                                 | Mains connection  | connector in acc. with EN 60320                         |  |
|                                 | Pump connection   | socket outlet in accordance with EN 60320               |  |
| Ambient Conditions              | Ambient temperature   | +10 to +40 °C for operation and storage                 |  |
|                                 | Air humidity  | max. 80 % in acc. with EN 61010, Part 1                 |  |
| Housing                         | Material  | steel aluminium housing                                 |  |
|                                 |   | with chemically resistant 2-component coating           |  |
|                                 | Dimensions  | (W x H x D) ~255 x 205 x 320 mm                         |  |
| Weight (incl. ViscoPump module) | ~5.4 kg   |   |  |
| Power supply                    | 90 to 240 V ~, 50 to 60 Hz  |   |  |
| Equipment safety                | EMC in acc. with Council Directive 89/336/EWG;<br>low-voltage directive |   |  |

## Ordering information AVS® 470

The AVS® 470 viscosity test station is composed of individual components. Please request a detailed quote..

| Type no.  | Order no. | Description  |
|---|-----------|--|
| AVS® 470 basic unit for opto-electronic sensing | 285415709 | AVS® 470 basic unit, housing incl. one ViscoPump III module for opto-electronic detection, keyboard<br>Version: 95 V to 230 V/50-60 Hz |
| AVS® 470 basic unit for TC sensing              | 285415708 | AVS® 470 basic unit, housing incl. one ViscoPump III module for TC detection, keyboard<br>Version: 95 V bis 230 V/50-60 Hz             |
| VZ 8561   | 285424060 | ViscoPump III module for optical detection   |
| VZ 8562   | 285424070 | ViscoPump III module for TC detection  |
| Z 910   | 285225640 | RS232-C data printer   |

The AVS® 470 allows the use of the following SI Analytics®-viscometers:  
Ubbelohde viscometer to DIN, Ubbelohde viscometer to ASTM, micro Ubbelohde viscometer to DIN, micro Ostwald viscometer, Cannon-Fenske routine viscometer, TC Ubbelohde viscometer, TC micro Ubbelohde viscometer.

We reserve the right to make technical changes.

AVS® is a registered trademark of SI Analytics® and stands for: "Automatic Viscosity System".

# 1.5 AVS® 370

## AVS® 370 makes maximum precision ...

### Well equipped for all viscosity determination

The AVS® 370 is a PC-controlled measuring device, which not only measures as precisely and consistently as you expect, but also offers maximum flexibility and future extensions. Furthermore, it saves laboratory space.

#### ▲ Suction and pressure mode - with one device

The AVS® 370 is operating with the ViscoPump III as the control unit for measurement and rinsing and therefore is able to pump the sample liquid in two different ways: by "suction" as well as by "pressure".

This makes it possible to simple adapt the method of measurement to different samples and applications: E.g. for non-critical samples as pharmaceutical solutions, we recommend to use the classic pressure mode. In pressure mode, at first the viscometer is filled and afterwards the filling and venting tube of



## ... easier and more flexible, with provision for future expansion!

the viscometer are connected to the ViscoPump. By applying pressure, the liquid is pumped to the measuring bulb.

The pressure mode is recommended for high-volatile solvents, as evaporation is lower compared to solvent mode. The viscometers have to be discharged and cleaned manually.

In suction mode, the sample liquid is pumped by vacuum up through the capillary. One main advantage is, that liquid can not leave the system. Therefore this method is often applied for hazardous samples.

For such samples, e.g. in polymer analytics, it is often requested to reduce liquid handling as much as possible, and therefore also the manual cleaning of viscometers. For these applications we recommend a waste system, combined with the AVS® 370, which makes manual discharge and cleaning of viscometers obsolete.

When combined with a waste system, the AVS® 370 exclusively is working in suction mode: due to increased safety, and the fact that in suction mode the filling tube of viscometer keeps open, which is required for sample filling.



- ▶ Automatic and highly precise measurements
- ▶ "Suction" and "pressure" measurements with the same module
- ▶ Modular concept for up to four ViscoPump III modules in one AVS® 370
- ▶ Each ViscoPump III module in a AVS® 370 can measure a different sample using a different method.
- ▶ Real multitasking for up to eight parallel measurements with the software Win-Visco 4
- ▶ TC version for measurement of nontransparent and black liquids

Advantages  
AVS® 370

As a further advantage, the suction mode exhibits increased reproducibility of flow times in case of Ubbelohde viscometers and samples which tend towards foaming: The bubbles are created during pumping of the liquid to the measuring bulb, when air is mixed with the small amount of liquid remaining in the lowest part of the capillary. In suction mode, the liquid can be blown out by a special function of software WinVisco 4.

## Two detection methods to measure flow times

To measure the flow time, the liquid meniscus can be detected by optoelectronics or thermally, by TC sensors. In both cases, the flow time is displayed with an resolution of 0.01 s.

Using optoelectronic detection, the liquid meniscus is registered by using IR light barriers; for detection with TC sensors, the different thermal conductivity of air and sample is utilized. With these options, the AVS370 covers a broad range of applications, including transparent and opaque samples as well.

### Easy modular concept ideal for future expansion

The AVS® 370 has a modular design. The basic version is available with one ViscoPump III module in optical or in TC version. Up to 3 other ViscoPump III modules can be installed in the compact housing. The measuring station can be adapted to increasing requirements at any time.

### Can be expanded from an affordable single measuring station up to an 8-sample station

The basic version of the AVS® 370 is able to measure the viscosity of liquids automatically. The TC version viscometers, it is ideal for measuring opaque and

black fluids. If necessary, each single measuring station can be expanded to a multiple measuring station with up to eight measuring positions. The WinVisco 4 software included with the standard equipment enables parallel operation of two fully equipped AVS® 370, with a total of eight ViscoPump III modules. Each module can measure a different sample using its own method. All the results can be quickly and easily evaluated and documented independently. It could hardly be more flexible!

### Waste system and rinsing

As mentioned above for the suction mode operation, the AVS® 370 can be combined with a waste system.

When using a waste system, after measurement the sample is discharged from the viscometer into a waste bottle - the viscometer is rinsed while keeping installed. The manual cleaning of the viscometer becomes obsolete, and the effort for dis- and reassembling to the measuring stand as well.



For 1 and 2 measuring positions, we offer the AVS®370 as a package, containing all components including measuring stands, thermostat bath, recirculating cooler, safety sensors, filling and waste system with discharge pump, complete hose sets and PC software.

▶ Compatible with existing accessories

Possibly existing accessories (thermostats, stands, flow through cooler, etc. von Vorgängermodellen) can continue to be used with the AVS® 370. Also, virtually all customary SI Analytics®-capillary viscometers can be used.



## AVS® 370 – the right solution for all situations

### Working with AVS® 370 is easy

The entire measurement procedure is placed automatic, subjective measuring errors are reliably eliminated. The PC starts the measurement. After the set preconditioning period the selected number of Durchflusszeiten gemessen and the measured values saved.

The system protects against accidental overpumping or oversuction by means of a capacitive sen-

sor. This prevents the sample to be measured from getting into the vessel containing the liquid or inside the device.

### Unique flexibility

In the PC-controlled multiple measurement station, the AVS® 370 offers unique flexibility while working in a very small space: Up to eight ViscoPump-modules, which equates to two fully equipped AVS® 370, can be run in parallel with the WinVisco 470 software.

Each module can measure the same or different samples using "pressure" or "suction", independent of each other. This significantly reduces the time required to carry out viscosity measurements in Parallelbetrieb, especially for in process controls and quality assurance. In this way, a series of measurements can be prepared quickly and immediately evaluated and documented with the computer.



## Technical data

|                                      |   |  |
|--------------------------------------|---|--|
| Measuring range (time)               | up to 9,999.99 s; resolution 0.01 s   |  |
| Measuring range (viscosity)          | pressure:   | 0.35 to 1,800 mm <sup>2</sup> /s (cSt)                                       |
|                                      | suction:  | 0.35 to ~5,000 mm <sup>2</sup> /s (cSt)                                      |
| Measured parameter                   | flow through time [s]   |  |
| Accuracy of the time measurement     | ± 0.01 %  |  |
| Measured value display               | via PC  |  |
| Display accuracy                     | ± 1 digit (0.1 %)   |  |
| Pump pressure                        | automatically controlled  |  |
| Preselectable tempering period       | 0 to 20 min   |  |
| Preselectable number of measurements | up to 10  |  |
| Connections                          | Pneumatic connections   | threaded connections for viscometers   |
|                                      | Electrical connections  | circular connector with bayonet lock for measuring stands and TC viscometers |
|                                      | RS232-C interface   | 9-pin  |
|                                      | Mains connections   | plug in accordance with EN 60320   |
|                                      | Pump connection   | socket outlet in accordance with EN 60320                                    |
| Data Input/Output                    | serial to EIA RS232-C   |  |
| Ambient conditions                   | Ambient temperature   | + 10 to + 40 °C  |
|                                      | Air humidity  | max. 85 % rel.   |
| Housing                              | Material  | coated aluminum plate  |
|                                      | Dimensions (for 1 to 4 modules)   | (W x H x D) ~255 x 205 x 320 mm  |
|                                      | Weight (incl. 1 module)   | ~5.4 kg  |
| Power supply                         | 90 to 240 V ~, 50 to 60 Hz  |  |
| Equipment safety                     | EMC-Compatibility according to the Directive 89/336/EEC of the Council<br>low-voltage directive according to the Directive 73/23/EEC of the Council<br>as amended by the Directive 93/68/EEC of the Council |  |
| Multi-tasking                        | for 1 to 8 ViscoPump III modules, with WinVisco 4 software  |  |

## Ordering information AVS® 370

The AVS® 370 viscosity test station is composed of individual components.  
Please request a detailed quote.

| Typ-Nr.   | Bestell-Nr. | Beschreibung  |
|---|-------------|---|
| AVS® 370 basic unit for opto-electronic detection | 1056509     | AVS® 370 basic unit, housing incl. one ViscoPump III-Moduls and Software WinVisco 370 software, for opto-electronic sensing |
| AVS® 370 basic unit for TC detection              | 1056515     | AVS® 370 basic unit, housing incl. one ViscoPump III-Moduls and Software WinVisco 370 software, for TC sensing              |

Accessories for AVS 370® and AVS 470® you can find at page 49 and following.

The following viscometers can be used with the AVS® 370:  
Ubbelohde viscometer to DIN, Ubbelohde viscometer to ASTM, micro Ubbelohde viscometer to DIN, micro Ostwald viscometer, Cannon-Fenske routine viscometer, TC-Ubbelohde viscometer, TC-micro Ubbelohde viscometer.

We reserve the right to make technical changes.

AVS® is a registered trademark of SI Analytics and stands for: "Automatic Viscosity System".

# 1.6 Software WinVisco 4

## Software WinVisco 4 ...



Clearly separated: work area, main menu, menu bar and header

The new WinVisco 4 is the ideal software to control the AVS®370, and it is included in the basic configuration of the instrument. Up to 8 viscometry measuring units (ViscoPumps) can be controlled by only few operating steps.

The WinVisco 4 is newly programmed - on the one hand combining clear representation of results, while on the other hand a comprehensive and detailed configuration of settings for the measuring method is possible.

As the previous version, the software is working in real multi-tasking operation, as WinVisco 4 is cooperating with the Visco-Pumps, whose internal software is controlling the time measurement and the measuring sequences.

### Individual configuration of user interface

- Multi-Language: English, German, French, Spanish, Chinese
- Changeable Layout, e.g. graphics in white/black or black/silver.



Individual Layout

## ..the new control software for AVS®370

### Simple operation

The daily measuring routine is carried out in a simple manner, using two windows: "Start" and "Overview"

### Temperature monitoring

New: With WinVisco 4, the bath temperature can be monitored when using a thermostat of the CT52 and CT72 series.

### Clearly arranged

The requested results, e.g. the kinematic viscosity or viscosity number, can be shown in the overview.

### User management

The WinVisco 4 supports three different user types of users. In the lowest level, user, only operation is allowed. It includes the execution of measurements, but also the selection of methods and viscometers. In addition, the lab manager can change all settings: Configuration of the measuring position and method as well as the viscometer database.

The administrator additionally can manage the user access data.

Sample input in window "Start"




Everything under control


| No. | Time [s] |
|-----|----------|
| 1   | 91.16    |
| 2   | 91.16    |

| Name                    | Value |
|-------------------------|-------|
| Rel. Viscosity          | 1.72  |
| Viscosity number [ml/g] | 144   |


Overview of current measurements



Configuration



Method management



User-defined formula

## Configuration of measurement station

For individual configuration of single measuring positions, e.g. "suction" or "pressure" mode, or activation of safety sensors.

## Method

All measurement settings which valid for all positions are done in the method, e.g. parameters of the ViscoPump III and the measurement, the evaluation and the configuration of an optional rinsing sequence.

## Rinsing

To rinse, there are two possibilities: rinsing with sample and rinsing with solvent. In both cases, the AVS® 370 is connected to a waste system, to enable rinsing of the viscometers keeping installed: With a built-in vacuum pump, samples and rinsing solvent are sucked into a waste bottle - the disassembling of viscometers for cleaning becomes obsolete.

- ▶ Simple handling in routing operation
- ▶ Individually configurable
- ▶ Multi-Language
- ▶ Monitoring the bath temperature with thermostat CT72
- ▶ Control up to eight measuring positions

Benefits  
WinVisco 4

The rinsing with sample is especially used in case of some polymer applications. Alternatively, rinsing with up to two solvents for each measuring position is possible.


In case of solvent rinsing, a daisy chain connection enables the integration and control of additional burettes Titronic® 300 (or Titronic® 500) to the AVS® 370.

▶ The sequences for rinsing are individually configured in the software.


## Evaluation

For evaluation, the WinVisco 4 includes many predefined formulas at state of delivery.

Additional calculations can be defined using a formel editor.



Clearly structured: Rinsing sequence



Transparent thermostat with four measuring stands

# 1.7 AVS® Pro III

## Automatic viscosity measurement has been improved ...

The AVS® Pro III automatic sampler is a fully automated instrument for determining the viscosity of Newtonian fluids using capillary viscometers. The system is mainly used for polymer analytics, e.g. according to ISO 307, ISO 1628 or ASTM D4603. But it also complies to the requirements of absolute measurements according to DIN 53 000, ASTM D445/D446 and ISO 3104/3105. In spite of the high sample throughput, the AVS® Pro III provides maximum accuracy and reproducibility. Furthermore, working with the automatic sampler is easy and allows unattended 24-hour operation.

In comparison to AVS® 370, at AVS® Pro III also the filling of the viscometer is automated. Therefore the complete measuring sequence - filling, measuring, discharging and rinsing - is fully automated, to enable unattended operation. The filling is done by a dosing module (piston/cylinder) in a 4-axis robotic system:

- ▶ Fully automatic and highly precise measuring station. Time measurement with a precision of +/- 0.01 s (but less precise than 0.1%)
- ▶ Ideal for highly aggressive media
- ▶ Although in combination of optical and thermal sampling of the meniscus channel or different capillary sizes and types, up to four viscometers selectable
- ▶ Durch die optionale Filtration mit dem ProClean system eliminate manual filtration of the sample
- ▶ Individually configurable, e.g. for rinsing with sample or solvent, or use of different viscometer types.

Advantages  
AVS® Pro III

The samples are sucked by a sample needle from a sample bottle and transferred to the viscometer. Due to medium hose lines of only minimum length, any carryover can be kept low easily.

Particularly for high sample throughput, the AVS® Pro III helps to substantially reduce the burden on qualified users. The fully automatic mode offers an increased level of safety when handling aggressive media, e.g. sulphuric acid.

### Applications

For the concept of AVS® Pro III, particularly the requirements of polymer quality control have been incorporated, e.g. for measurement of intrinsic viscosity, viscosity number, K value, or similar measures. But also the measurement of petrochemical products is possible.

### Different configurations

The AVS® Pro III typically is equipped with 2 to 4 measuring positions, to achieve higher sample throughput when measuring in parallel.

The following viscometer types can be used: Ubbelohde (DIN, ASTM), Mikro-Ubbelohde (DIN), Cannon-Fenske routine and Micro-Ostwald of SI Analytics®.

The standard waste system contains a waste bottle for each measuring position. The advantage: This system is very reliable in operation, as it does not need valves being in contact to the medium. But as a special configuration, also a waste system with only one waste bottle is available.

## ... with the AVS® Pro III Automatic Sampler

### Safety and reliability

The AVS® Pro III has an outstanding reliability: The dosing system (available as standard and Micro version) is operating without valve and therefore is suitable for nearly all samples. Due to further development of hardware components and software, the system has been continuously improved, resulting in a minimization of error sources. E.g., in the robotic system, proximity switches are in use. And the system abstains from mediumconveying valves, if not required for special customer applications. Even more

important than reliability is the safety for the user. For this, the liquid in the waste system is transferred only by suction instead of high pressure: Therefore, the liquid is only aspirated to the waste bottle, not pressed. In case of leakage, only air can enter into the system, but never liquid leave the system.

Regarding safety devices, there are capacitive sensors to control the suction line between viscometer and control unit, as well as liquid level detection of the waste bottles. In case of a signal of these safety devices, the corresponding measuring position is switched off.



## Rinsing: Either with sample or with solvent

The rinsing of dosing module and viscometer is carried out according to the application and customer's preferences. Both alternatives have benefits and drawbacks.

Most AVS® III systems are configured for rinsing with next sample. For this, in total 3 fillings of viscometer are required: Using the first 2 charges, the dosing module and the viscometer are rinsed, and only the 3rd filling is used for measurement.

As Ubbelohde viscometers of standard size require a liquid volume of about 17 ml, the total sample volume needed is about 60 ml.


Therefore typically sample bottles of 100 ml volume are used. We recommend this configuration for applications in polymer analytics, when supply and disposal of the sample solvent is cheap and rinsing with external solvent exhibits drawbacks. An important application is the measurement of polyamide in sulfuric acid (96%) according to ISO 307.

When using expensive and hazardous solvents, e.g. for measurement of polyesters according to ISO 1628-5 or ASTM D4603, often the user wants to use only a small amount of sample. In this case it is possible to use micro viscometers with filling volume of max. 4 ml - in this case a sample volume of 20 ml is sufficient for rinsing and measurement. In case of micro viscometers, the AVS® Pro III is equipped with a special dosing module for small sample volumes.

Alternatively, instead changing to viscometers with smaller sample volume, there is the option of rinsing with external solvent - then only one filling of viscometer is required. In this way, also for standard Ubbelohde viscometers a total sample volume of 25 ml is sufficient. Using solvent rinsing, dosing module and viscometer are rinsed 2 times with solvent. Preferably, this solvent is volatile, to enable a drying step directly afterwards. Otherwise, subsequent a second solvent with low boiling point has to be used, but this does not need to be a solvent for original sample. During drying the remaining residues of

solvent are evaporated by applying vacuum respective in a stream of air, therefore afterwards the next sample can be filled in.

When rinsing with sample, not only some reagent for sample preparation is saved: due to smaller sample bottles, sample racks with 56 instead of 16 positions can be used.



4-tube viscometer, for rinsing with solvent

As a special option it is also possible, to insert an additional rinsing sequence at the end of a measuring series.

### Sample rack


There are several sample racks available, depending on the size of the sample bottles. For 100 ml bottles a sample rack with 16 positions is used. For 40 ml bottles and 20 ml vials, there is a sample rack with 56 positions available. The samples are positioned in the sample rack, which is easy to load due to an electric sample lift. The AVS® Pro III allows free selection of sample sequence and moreover individual allocation of samples to certain viscometers.



Two different sample racks are supplied:  
a) rack with 56 positions for 20 ml and 40 ml (new) sample bottles.




b) Rack with 16 positions for 100ml sample bottles and viscometer with standard sample volume and rinsing with sample



The electric sample lift ensures positioning of the samples in the rack at a convenient and easily monitored working height.

## Filtration

As long as particle-free liquids are measured, there is no filtration necessary. In case of only occasional particles, inline filters can be integrated into the filling hoses, to avoid any issues by impurities entering the viscometer. But in case of measuring e.g. glass fiber infoced polymer samples, inline filter can't be used - they would be blocked within short time of use. For such samples there is the ProClean filtration system available, which is already filtering the sample in the sample bottle.



Filtration system ProClean

## Working with the AVS® Pro III is ...

The AVS® Pro III is controlled by a PC. The intuitive user interface guides the user clearly through the program. All data inputs are made using the computer keyboard and mouse.

A faulty operating status is indicated by acoustic or optical signals such as arrows, icons and other status messages or request messages. During the entire work sequence, the respective status of the AVS® Pro III is documented on the computer screen. Furthermore, status indicators can be selected for each individual measuring position, which provide additional information on operation.

For the respective type of measurement, pre-parameterized sets of parameters depending on the viscometers, temperature and other measurement criteria are already provided. In addition, all parameters can be individually adjusted to special requirements at a special menu level. All of the standard calculation methods are available:

- mean value
- standard deviation
- outlier test (A %)
- Hagenbach correction
- absolute viscosity, dynamic viscosity (density value required)
- viscosity index (measurement at two temperatures required)
- SUS and SFS
- relative viscosity
- specific viscosity
- reduced viscosity (viscosity number)
- inherent viscosity
- intrinsic viscosity and
- K-value after Fikentscher

The proved and tested AVS® Pro III software also makes it possible to prepare additional individually selected calculations:

During the entire process, all of the parameters (depending on the menu level) and the respective status of the individual measuring positions, the temperature regulation system and the sample transfer system are either visible or can be selected.

The operator interface of the AVS® Pro III is available in German and English. Commercially available

printers for which Windows drivers are available are suitable for documentation purposes.

The AVS® Pro III is built in accordance with international equipment safety standards: and CE certified (equipment safety, low voltage safety, emitted interference and interference immunity).


If requested, the AVS® Pro III automatic sampler can be supplied with a manufacturer's inspection certificate based on direct comparison with normal viscometers of the first order in accordance with DIN 53 000-3.

### Software

Some features of the software are described in the following.

#### Individual sample allocation

The AVS® Pro III allows individual allocation between the characteristics of the sample and the viscometers that are currently in operation.



*The allocation between the sample and the viscometer is shown on the status display.*

In practice, samples of different viscosities can be filled into viscometers of different size and measured at the same time. This even applies to a combination of optical and thermal scanning. Therefore, preliminary sorting of the samples with regard to viscosity and the size of capillary required for the testing process is not necessary.

## ... easy, reliable and safe

### Method

This mode is used to specify what monitoring parameters are to be activated, e.g. if the temperature control of the thermostats is supposed to be handled via the PC.

### Options

In several windows, configuration details are parametrized: e.g. size of dosing module, immersion depth of the sample needle, volumes of rinsing solvent, activation of safety devices and temperature control of bath thermostat, but also change of language English/German etc.

### Selection of method

| Immersion depth  |        |
|------------------|--------|
| Sample           | 123 mm |
| solvent 1        | 22 mm  |
| solvent 2        | 22 mm  |
| Waste port       | 25 mm  |
| Viscometer ports | 30 mm  |

### Options



# Technical data AVS® Pro III

|                               |   |   |  |
|-------------------------------|---|---|--|
| Sampling system               | Sample bottles  | 100 ml screw-type and bottles with standard ground joint (16 pcs per rack)          |  |
|                               |   | 20 ml round bottom glass pieces (56 pcs. per rack)                                  |  |
|                               |   | 40 ml EPA-Screw-thread bottle (16 or 56 pcs. per rack)                              |  |
|                               | Sample rack   | for 100 ml screw-type and bottles with standard ground joint                        |  |
|                               |   | for 20 ml round bottom glass pieces   |  |
|                               |   | for 40 ml EPA-Screw-thread bottle (16 or 56 pcs. per rack)                          |  |
| Measured value recording      | Method  | meniscus scanning by means of opto-electronic system or thermal conductivity (TC)   |  |
| Measuring parameter           | throughput time in seconds [s]  |   |  |
|                               | temperature in degrees Celsius [°C]   |   |  |
| Calculated parameters         | mean value, standard deviation, outlier test (A %), Hagenbach correction, absolute viscosity, dynamic viscosity (knowledge of density required), viscosity index (measurement at two temperatures required) SUS and SFS, relative viscosity, specific viscosity, reduced viscosity (viscosity number), inherent viscosity, K-value, intrinsic viscosity   |   |  |
| Selection parameters          | by means of PC keyboard, mean value, standard deviation, outlier test (A %), Hagenbach correction, absolute viscosity, dynamic viscosity (knowledge of density required), viscosity index (measurement at two temperatures required) SUS and SFS, relative viscosity, specific viscosity, reduced viscosity (viscosity number), inherent viscosity, K-value, rack position, date/time, temperature regulation period, number of measurements, number of rinsing operations, start, stop/reset |   |  |
|                               | Number of measurements  | 1 to 99   |  |
|                               | Temperature regulation period   | 0 to 99 min., selectable in increments of 1 min.                                    |  |
|                               | Number of Viscometer tests  | 0 to 10 with next sample (observe sample quantity) or with external rinsing solvent |  |
|                               | Data memory   | by means of PC  |  |
| Viscosity measurement range   | 0.35 to 1,200 mm <sup>2</sup> /s (at room temperature of samples)   |   |  |
|                               | Time  | up to 9999.99 s, resolution = 0.01 s  |  |
|                               | Vacuum pressure   | automatically controlled  |  |
|                               | Viscometers available for use   | Ubbelohde viscometer in accordance with DIN standards                               |  |
|                               |   | Ubbelohde viscometer in accordance with ASTM standards                              |  |
|                               |   | Micro-Ubbelohde viscometer in accordance with DIN standards                         |  |
|                               |   | Micro-Ostwald viscometer, Cannon-Fenske-Routine visco                               |  |
|                               |   | Cannon-Fenske-Routine viscometer  |  |
| TC Ubbelohde viscometer       |   |   |  |
| TC Micro-Ubbelohde viscometer |   |   |  |

|                       |   |  |
|-----------------------|---|--|
| Measuring accuracy    | ± 0.01 s ± 1 digit, but not more precise than 0.01 %  |  |
|                       | The measuring uncertainty for measurements of absolute kinematic viscosity is also dependent on the uncertainty of the numeric value for the viscometer constant and on the measuring conditions, especially the measuring temperature. |  |
| Evaluations / results | Correction  | Hagenbach correction (HC) for Ubbelohde, Cannon-Fenske-Routine, Micro-Ubbelohde and Micro-Ostwald viscometers      |
|                       | Statistical evaluation  | standard deviation, outlier search   |
| Ambient conditions    | Ambient temperature   | 10 to + 40 °C  |
|                       | Air humidity  | max. 85 % relative humidity  |
| Equipment safety      | CE-symbol   | in accordance with Guideline 89/336/EEC of the Council (EMC compatibility)   |
|                       |   | in accordance with Standard EN 50 081, Part 1  |
|                       |   | interference immunity in accordance with Standard EN 50 082, Part 2  |
|                       |   | in accordance with Guideline 73/23/EEC of the Council (low-voltage guideline)                                      |
| Housing               | plastic/stainless steel / aluminium casing with chemically resistant two-component coating of the plastic pieces  |  |
|                       | Dimensions  | w = 1.300 mm, h = 1.000 mm, d = 620 mm (~51" x 43" x 24")  |
|                       | Weight  | dependent on the number of measuring positions<br>~70 kg   |
| Connections           | Pneumatic connections   | screw-type connections for viscometer  |
|                       | Electric connections  | circular connectors with bayonet lock for measuring stand and TC viscometer  |
|                       | Viscometers   | up to 4 viscometers connected by individual control units  |
|                       | Temperature   | via serial interface RS232-C of suspended thermostat<br>type: 1 pc. CT 72/4 or up to 2 pcs. CT 72/2, 1 pc. CT 72/2 |
|                       | Interfaces  | control system using PC with 2 x RS232-C interfaces bzw. USB/RS232-Adapter   |
|                       | Safety  | overflowing safety device of waste bottle and suction hose   |
|                       | Mains connection  | European built-in plug DIN 49 457 6 with fuse  |
| Data transmission     | Interface internal  | bidirectional serial interface in accordance with EIA RS232-C (daisy chain concept)                                |
|                       | Interface external  | via PC, bidirectional serial interface in accordance with EIA RS232-C  |
| Power supply          | Mains voltage   | 230 V (AC) or 115 V (AC), 50 to 60 Hz (AC)   |

# 1.8 CT 72 Thermostat Series



CT 72/4

- ▶ CT 72/2 and CT 72/4 can be used up to 150 °C. High temperature version is standard.
- ▶ Draining valve comes with CT 72/2, CT 72/2-TT and CT 72/4.
- ▶ High temperature stability  $\pm 0.02$  K as requested in viscometry standards

Advantages  
Baths

As their predecessor CT 52 the transparent thermostats CT 72/P, CT 72/2, CT 72/2-TT and CT 72/4 meet DIN 53 000 part 1, ASTM D 445 and ISO 3105 standards.



CT 72/2

The SI Analytics® transparent thermostats are particularly designed for the determination of the viscosity of newtonian liquids in glass capillary viscometers. They may be adapted for manual as well as for automatic measurements.



Bright display shows ongoing process, at any time.

- ▶ Programmable set temperatures through integrated clock with controller.
- ▶ Display of the momentary and the set temperature.
- ▶ Increased safety from separate operation and temperature safety sensors.
- ▶ Automatic fuses on the back panel.

New  
Immersion  
thermostats

# Transparent Thermostats conforming to Standards: The CT 72 series



CT 72/P



CT 72/4

## Suitable temperature control liquids

| Liquid            | Alcohol          | Water           | Paraffine oil     | Silicon oil       |
|-------------------|------------------|-----------------|-------------------|-------------------|
| Temperature range | -40 °C to +10 °C | +5 °C to +80 °C | +40 °C to +150 °C | +80 °C to +150 °C |


## Technical data

| Device  | CT 72/P          | CT 72/2-TT        | CT 72/2           | CT 72/4           |
|---|------------------|-------------------|-------------------|-------------------|
| Working temperature                                 | +10 °C to +60 °C | -40 °C to +150 °C | +5 °C to +150 °C  | +5 °C to +150 °C  |
| Measuring positions for AVS                         | 2                | 2                 | 2                 | 4                 |
| Measuring positions TC                              | 2                | 2                 | 2                 | 4                 |
| Measuring positions micro TC                        | 2                | 2                 | 2                 | 4                 |
| Temperature stability according DIN 58 966 at 25 °C | ±0.02 K          | ±0.02 K           | ±0.02 K           | ±0.02 K           |
| Size (W x H x D in mm)                              | 355 x 370 x 250  | 355 x 370 x 250   | 355 x 370 x 250   | 605 x 370 x 250   |
| Filling volume                                      | 18 l             | 15 l              | 15 l              | 27 l              |
| Material  | PMMA             | St. steel & glass | St. steel & glass | St. steel & glass |
| Weight (empty)                                      | ~5 kg            | ~14 kg            | ~13.5 kg          | ~28 kg            |

At applications within normal temperature range (+5 °C up to approx. +40 °C) cooling will be necessary for maintaining temperature stability. This cooling can be carried out by flowing tap water or by use of an optional flow-through cooler (e.g. CK 310). For low temperature applications, an optional cryostat with high cooling power is required.

## Ordering information

| Type no.                         | Order no. | Description  |
|----------------------------------|-----------|--|
| CT 72/P, 230V                    | 285418526 | Immersion thermostat 230 V and thermostatic bath (acrylic glass container with two manual gauge slides), basic configuration for the attachment of one flow-through cooler.  |
| CT 72/P, 115V                    | 285418513 | Immersion thermostat 115 V and thermostatic bath (acrylic glass container with two manual gauge slides), basic configuration for the attachment of one flow-through cooler.  |
| CT 72/2, 230V                    | 285418547 | Immersion thermostat 230 V and thermostatic bath (stainless steel container with one manual gauge slide), basic configuration for the attachment of one flow-through cooler.   |
| CT 72/2, 115V                    | 285418532 | Immersion thermostat 115 V and thermostatic bath (stainless steel container with one manual gauge slide), basic configuration for the attachment of one flow-through cooler.   |
| CT 72/2-M, 230V                  | 285418584 | Immersion thermostat 230 V and thermostatic bath (stainless steel container with one manual gauge slide), equipped with two magnetic stirrer positions. Basic configuration for the attachment of one flow-through cooler. |
| CT 72/2-M, 115V                  | 285418593 | Immersion thermostat 115 V and thermostatic bath (stainless steel container with one manual gauge slide), equipped with two magnetic stirrer positions. Basic configuration for the attachment of one flow-through cooler. |
| CT 72/2-TT, 230V                 | 285418615 | Immersion thermostat 230 V and thermostatic bath (stainless steel container with one manual gauge slide), basic configuration for the attachment of one flow-through cooler.   |
| CT 72/2-TT, 115V                 | 285418607 | Immersion thermostat 115 V and thermostatic bath (stainless steel container with one manual gauge slide), basic configuration for the attachment of one flow-through cooler.   |
| CT 72/4, 230V                    | 285418568 | Immersion thermostat 230 V and thermostatic bath (stainless steel container with two manual gauge slides), basic configuration for the attachment of one flow-through cooler.  |
| CT 72/4, 115V                    | 285418554 | Immersion thermostat 115 V and thermostatic bath (stainless steel container with two manual gauge slides), basic configuration for the attachment of one flow-through cooler.  |
| CT 72/E, 230V                    | 285418501 | Immersion thermostat 230 V/50 Hz   |
| CT 72/E, 115V                    | 285418495 | Immersion thermostat 115 V/60 Hz   |
| CK 310, 230V                     | 285414320 | Flow-through cooler CK 310, version: 230 V/50-60 Hz with cooling coil made of stainless steel  |
| CK 310, 115V                     | 285414310 | Flow-through cooler CK 310, version: 115 V/50-60 Hz with cooling coil made of stainless steel  |
| More Accessories and spare parts |           |  |
| VZ 5402                          | 285415171 | Manual gauge slide for transparent thermostats   |
| VZ 5403                          | 285420684 | 3-fold manual gauge slide for transparent thermostats  |
| VZ 5404                          | 285418573 | Dust protection cover for transparent thermostat   |
| VZ 5405                          | 285418620 | Transparent thermostatic bath backlight  |



CT 72/2



CK 310

## 2 Viscometers

### 2.1 Viscometers and their range of use

| Measurement substance property                                | Viscometer type |                 |                 |         |               |                       |                            |
|---|-----------------|-----------------|-----------------|---------|---------------|-----------------------|----------------------------|
|   | Ubbelohde       | Micro-Ubbelohde | TC Ubbelohde    | Ostwald | Micro-Ostwald | Cannon-Fenske-Routine | Cannon-Fenske-reverse flow |
| Transparent liquids manual measurement                        | ✓               | ✓               | ■               | ■       | ■             | ■                     | ■                          |
| Transparent liquids automatic measurement                     | ✓               | ✓               | ■               | ■       | ■             | ■                     | ■                          |
| Opaque liquids manual measurement                             | ■               | ■               | ■               | ■       | ■             | ■                     | ■                          |
| Opaque liquids automatic measurement                          | ■               | ■               | ✓ <sup>1)</sup> | ■       | ■             | ■                     | ■                          |
| Foaming liquids   | ■               | ■               | ■               | ■       | ■             | ■                     | ■                          |
| Liquid mixture with highly volatile components                | ■               | ■               | ■               | ■       | ■             | ■                     | ■                          |
| Minimum measurement substance and/or rinsing agent quantities | ■               | ✓               | ■               | ■       | ✓             | ■                     | ■                          |
| High-temperature or low-temperature measurements              | ■               | ■               | ■               | ■       | ■             | ■                     | ■                          |

Selection of glass capillary viscometers

|   |                 |
|---|-----------------|
| ✓ | use preferably  |
| ■ | highly suitable |
| ■ | less suitable   |
| ■ | unsuitable      |

<sup>1)</sup> up to 30,000 mm<sup>2</sup>/s  
<sup>2)</sup> above 30.000 mm<sup>2</sup>/s

## 2.2 Ubbelohde viscometers normal form (DIN)

Viscometers with suspended level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 53 350, Part 18. All viscometers are provided with ring marks. This ensures that viscometers for automatic mea-

surements can also be checked by means of manual measurements. The recommended minimum flowthrough time is 200s for absolute measurements of kinematic viscosity. For relative measurements (polymer analytics), the minimum flow time of 50s is allowed by ISO 1628-1 (depending on the capillary size).



| calibrated, with constant, for manual measurements |           | calibrated, with constant for manual measurements; automatic measurement with stand AVS®/SK-HV |           | $v = K \cdot t$<br>$K = \frac{v}{t}$<br>$t = \frac{v}{K}$ | Ubbelohde-Viskosimeter (DIN)          |                      |  |
|--|-----------|--|-----------|---|---------------------------------------|----------------------|--|
| Type No.   | Order No. | Type No.   | Order No. | Capillary No. acc. DIN                                    | Capillary $\varnothing \pm 0,01$ [mm] | Constant K (approx.) | Measuring range [mm <sup>2</sup> /s] (approx.) |
| 501 00   | 285400004 | -  | -         | 0   | 0.36                                  | 0.001                | 0.3 to 1                                       |
| 501 03   | 285400012 | -  | -         | 0c  | 0.47                                  | 0.003                | 0.5 to 3                                       |
| 501 01   | 285400029 | -  | -         | 0a  | 0.53                                  | 0.005                | 0.8 to 5                                       |
| 501 10   | 285400037 | -  | -         | I   | 0.63                                  | 0.01                 | 1.2 to 10                                      |
| 501 13   | 285400045 | -  | -         | Ic  | 0.84                                  | 0.03                 | 3 to 30  |
| 501 11   | 285400053 | -  | -         | Ia  | 0.95                                  | 0.05                 | 5 to 50  |
| 501 20   | 285400061 | -  | -         | II  | 1.13                                  | 0.1                  | 10 to 100                                      |
| 501 23   | 285400078 | -  | -         | IIc   | 1.50                                  | 0.3                  | 30 to 300                                      |
| 501 21   | 285400086 | -  | -         | IIa   | 1.69                                  | 0.5                  | 50 to 500                                      |
| 501 30   | 285400094 | -  | -         | III   | 2.01                                  | 1                    | 100 to 1,000                                   |
| 501 33   | 285400107 | -  | -         | IIIc  | 2.65                                  | 3                    | 300 to 3,000                                   |
| 501 31   | 285400115 | -  | -         | IIIa  | 3.00                                  | 5                    | 500 to 5,000                                   |
| 501 40   | 285400123 | -  | -         | IV  | 3.60                                  | 10                   | 1,000 to 10,000                                |
| -  | -         | 502 43   | 285400131 | IVc   | 4.70                                  | 30                   | 3,000 to 30,000                                |
| -  | -         | 502 41   | 285400148 | IVa   | 5.34                                  | 50                   | 6,000 to 30,000                                |
| -  | -         | 502 50   | 285400156 | -   | 6.30                                  | 100                  | above 10,000                                   |

| not calibrated, without constant; for determination of relative viscosity |           | calibrated, with constant for automatic measurements |           | $v = K \cdot t$<br>$K = \frac{v}{t}$<br>$t = \frac{v}{K}$ | $v = \frac{K}{t}$                     | kinematic viscosity in mm <sup>2</sup> /s = constant [mm <sup>2</sup> /s]<br>= flow-through time in s |  |
|---|-----------|--|-----------|---|---------------------------------------|---|--|
| Type No.  | Order No. | Type No.   | Order No. | Capillary No. acc. DIN                                    | Capillary $\varnothing \pm 0,01$ [mm] | Constant K (approx.)  | Measuring range [mm <sup>2</sup> /s] (approx.) |
| -   | -         | 532 00   | 285400164 | 0   | 0.36                                  | 0.001   | 0.3 to 1                                       |
| 530 03  | 285400304 | 532 03   | 285400201 | 0c  | 0.47                                  | 0.003   | 0.5 to 3                                       |
| 530 01  | 285400312 | 532 01   | 285400218 | 0a  | 0.53                                  | 0.005   | 0.8 to 5                                       |
| 530 10  | 285400329 | 532 10   | 285400226 | I   | 0.63                                  | 0.01  | 1.2 to 10                                      |
| 530 13  | 285400337 | 532 13   | 285400234 | Ic  | 0.84                                  | 0.03  | 3 to 30  |
| 530 11  | 285400338 | 532 11   | 285400172 | Ia  | 0.95                                  | 0.05  | 5 to 50  |
| 530 20  | 285400345 | 532 20   | 285400242 | II  | 1.13                                  | 0.1   | 10 to 100                                      |
| 530 23  | 285400353 | 532 23   | 285400259 | IIc   | 1.50                                  | 0.3   | 30 to 300                                      |
| 530 21  | 285400350 | 532 21   | 285400189 | IIa   | 1.69                                  | 0.5   | 50 to 500                                      |
| 530 30  | 285400361 | 532 30   | 285400267 | III   | 2.01                                  | 1   | 100 to 1,000                                   |
| 530 33  | 285400378 | 532 33   | 285400275 | IIIc  | 2.65                                  | 3   | 300 to 3,000                                   |
| 530 31  | 285400370 | 532 31   | 285400197 | IIIa  | 3.00                                  | 5   | 500 to 5,000                                   |
| 530 40  | 285400386 | 532 40   | 285400283 | IV  | 3.60                                  | 10  | 1,000 to 10,000                                |

## 2.3 Ubbelohde viscometers normal form (ASTM)



Viscometers with suspended level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18. All viscometers are provided with ring marks. This ensures that viscometers for automatic mea-

surements can also be checked by means of manual measurements. The recommended minimum flowthrough time is 200s for absolute measurements of kinematic viscosity. For relative measurements (polymer analytics), the minimum flow time of 50s is allowed by ISO 1628-1 (depending on the capillary size).

| calibrated,<br>with constant and<br>calibration<br>certificate for<br>manual<br>measurements |           | not calibrated,<br>without calibration<br>certificate for<br>determination of<br>relative Viscosity |           | calibrated,<br>with constant and<br>calibration<br>certificate for<br>automatic<br>measurements |           | Capillary<br>No. | Capillary<br>Ø $i \pm 0,01$<br>[mm] | Constant K<br>(approx.) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|--|-----------|---|-----------|---|-----------|------------------|-------------------------------------|-------------------------|--|
| Type<br>No.  | Order No. | Type<br>No.   | Order No. | Type<br>No.   | Order No. |                  |                                     |                         |  |
| 525 00   | 285400501 | 526 00  | 285400707 | 527 00  | 285401255 | 0                | 0.24                                | 0.001                   | 0.35 to 1  |
| 525 03   | 285400518 | 526 03  | 285400715 | 527 03  | 285401271 | 0C               | 0.36                                | 0.003                   | 0.6 to 3   |
| 525 01   | 285400526 | 526 01  | 285400723 | 527 01  | 285401263 | 0B               | 0.46                                | 0.005                   | 1 to 5   |
| 525 10   | 285400534 | 526 10  | 285400731 | 527 10  | 285401152 | 1                | 0.58                                | 0.01                    | 2 to 10  |
| 525 13   | 285400542 | 526 13  | 285400748 | 527 13  | 285401169 | 1C               | 0.78                                | 0.03                    | 6 to 30  |
| 525 11   | 285400550 | 526 11  | 285400750 | 527 11  | 285401170 | 1B               | 0.88                                | 0,05                    | 10 to 50   |
| 525 20   | 285400559 | 526 20  | 285400756 | 527 20  | 285401177 | 2                | 1.03                                | 0.1                     | 20 to 100  |
| 525 23   | 285400567 | 526 23  | 285400764 | 527 23  | 285401185 | 2C               | 1.36                                | 0.3                     | 60 to 300  |
| 525 30   | 285400575 | 526 30  | 285400772 | 527 30  | 285401193 | 3                | 1.83                                | 1                       | 200 to 1,000   |
| 525 33   | 285400583 | 526 33  | 285400789 | 527 33  | 285401288 | 3C               | 2.43                                | 3                       | 600 to 3,000   |
| 525 40   | 285400591 | 526 40  | 285400797 | 527 40  | 285401296 | 4                | 3.27                                | 10                      | 2,000 to 10,000                                      |
| 525 43   | 285400604 | 526 43  | 285400801 | 527 43  | 285401309 | 4C               | 4.32                                | 30                      | 6,000 to 30,000                                      |


### Ubbelohde Viscometer (ASTM)

- in accordance with ISO 3105, ASTM D 446
- filling quantity: 15 to 20 ml
- overall length: approx. 285 mm

## 2.4 Ubbelohde viscometers, with additional tube and threads

Viscometers with suspended level for determination of absolute or relative kinematic viscosity. These viscometers are preferably used for automatic measurements. The additional filling and cleaning tube and the glass thread ensure safe operational use. The cali-

brated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18. The ring marks present serve as auxiliary marks in case the viscometers must be checked by manual measurements.



calibrated,  
with constant and  
calibration certificate  
for automatic  
measurements

### Ubbelohde viscometer (DIN)

- in accordance with ISO 3105, DIN 53 000 Part 1
- filling quantity: 18 to 22 ml
- overall length: approx. 290 mm

| Type No. | Order No. | Capillary No. | Capillary Ø i [mm] | Constant K (approx.) | Measuring range [mm <sup>2</sup> /s] (approx.) |
|----------|-----------|---------------|--------------------|----------------------|--|
| 541 03   | 285401925 | 0c            | 0.47               | 0.003                | 0.5 to 3                                       |
| 541 01   | 285401917 | 0a            | 0.53               | 0.005                | 0.8 to 5                                       |
| 541 10   | 285401933 | I             | 0.63               | 0.01                 | 1.2 to 10                                      |
| 541 13   | 285401941 | Ic            | 0.84               | 0.03                 | 3 to 30  |
| 541 11   | 285401950 | Ia            | 0.95               | 0.05                 | 5 to 50  |
| 541 20   | 285401958 | II            | 1.13               | 0.1                  | 10 to 100                                      |
| 541 23   | 285401966 | IIc           | 1.50               | 0.3                  | 30 to 300                                      |
| 541 21   | 285408719 | IIa           | 1.69               | 0,5                  | 50 to 500                                      |
| 541 30   | 285401974 | III           | 2.01               | 1                    | 100 to 1,000                                   |
| 541 33   | 285401982 | IIIc          | 2.65               | 3                    | 300 to 3,000                                   |
| 541 40   | 285401999 | IV            | 3.60               | 10                   | 1,000 to 10,000                                |
| 541 43   | 285402000 | IVc           | 4,70               | 30                   | 3.000 to 30,000                                |


calibrated,  
with constant and calibration certificate  
for automatic measurements

### Ubbelohde viscometer (ASTM)

- the technical measurement characteristics are in accordance with ISO 3105, ASTM D 446
- filling quantity: 15 to 22 ml
- overall length: approx. 290 mm

| Type No. | Order No. | Capillary No. acc. ASTM | Capillary Ø i [mm] | Constant K (approx.) | Measuring range [mm <sup>2</sup> /s] (approx.) |
|----------|-----------|-------------------------|--------------------|----------------------|--|
| 545 00   | 285402005 | 0                       | 0.24               | 0.001                | 0.35 to 1                                      |
| 545 03   | 285402021 | 0C                      | 0.36               | 0.003                | 0.6 to 3                                       |
| 545 01   | 285402013 | 0B                      | 0.46               | 0.005                | 1 to 5   |
| 545 10   | 285402038 | 1                       | 0.58               | 0.01                 | 2 to 10  |
| 545 13   | 285402046 | 1C                      | 0.78               | 0.03                 | 6 to 30  |
| 545 11   | 285402042 | 1B                      | 0.88               | 0.05                 | 10 to 600                                      |
| 545 20   | 285402054 | 2                       | 1.03               | 0.1                  | 20 to 100                                      |
| 545 23   | 285402062 | 2C                      | 1.36               | 0.3                  | 60 to 300                                      |
| 545 30   | 285402079 | 3                       | 1.83               | 1                    | 200 to 1,000                                   |
| 545 33   | 285402087 | 3C                      | 2.43               | 3                    | 600 to 3,000                                   |
| 545 40   | 285402095 | 4                       | 3.27               | 10                   | 2,000 to 10,000                                |
| 545 43   | 285402108 | 4C                      | 4.32               | 30                   | 6,000 to 30,000                                |

## 2.5 Ubbelohde viscometers with TC sensors



Viscometers with suspended level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. The measuring levels are marked by TC sensors. The meniscus passage is detected due to the different conductivity of the liquid phase and the gas phase. A measurement stand of the type series AVS®/S is not required. TC viscometers can be used to determine the kinematic viscosity of all liquids with Newtonian flow behavior.


They are especially suitable for liquids that cannot be detected with other systems: opaque and/or black samples.

Due to the electrical properties of TC sensors, it is important to ensure that a suitable type is selected for the required application temperature.

### TC viscometers with additional filling and cleaning tube and with glass thread

- the technical measurement characteristics are in accordance with DIN 53 000, part 1, ISO 3105
- for use in combination with an automatic viscosity measuring instrument
- filling quantity: 18 to 22 ml
- overall length: approx. 355 mm
- suitable bracket Type No. 05393, Order No. 285405035

| calibrated,<br>with constant and calibration certificate for automatic<br>measurements |           |                 |           |                  |           | Capillary<br>No. | Capillary<br>Ø i [mm] | Constant K<br>(approx.) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|--|-----------|-----------------|-----------|------------------|-----------|------------------|-----------------------|-------------------------|--|
| Type<br>No.  | Order No. | Type<br>No.     | Order No. | Type<br>No.      | Order No. |                  |                       |                         |  |
| + 10 to + 80 °C  |           | - 40 to + 30 °C |           | + 70 to + 150 °C |           |                  |                       |                         |  |
| 562 03   | 285423120 | -               | -         | -                | -         | 0c               | 0,47                  | 0,003                   | 0,5 to 3   |
| 562 10   | 285423130 | 563 10          | 285423240 | 564 10           | 285423330 | I                | 0,63                  | 0,01                    | 1,2 to 10  |
| 562 13   | 285423140 | 563 13          | 285423250 | 564 13           | 285423340 | Ic               | 0,84                  | 0,03                    | 3 to 30  |
| 562 20   | 285423150 | 563 20          | 285423260 | 564 20           | 285423350 | II               | 1,13                  | 0,1                     | 10 to 100  |
| 562 23   | 285423170 | 563 23          | 285423270 | 564 23           | 285423360 | IIc              | 1,51                  | 0,3                     | 30 to 300  |
| 562 21   | 285423160 | -               | -         | -                | -         | IIa              | 1,69                  | 0,5                     | 50 to 500  |
| 562 30   | 285423180 | 563 30          | 285423280 | 564 30           | 285423370 | III              | 2,05                  | 1                       | 100 to 1.000   |
| 562 33   | 285423200 | 563 33          | 285423290 | 564 33           | 285423380 | IIIc             | 2,7                   | 3                       | 300 to 3.000   |
| 562 31   | 285423190 | -               | -         | -                | -         | IIIa             | 3,0                   | 5                       | 500 to 5.000   |
| 562 40   | 285423210 | 563 40          | 285423300 | 564 40           | 285423390 | IV               | 3,7                   | 10                      | 1.000 to 10.000                                      |
| 562 43   | 285423230 | 563 43          | 285423320 | 564 43           | 285423400 | IVc              | 4,9                   | 30                      | 3.000 to 20.000                                      |
| 562 41   | 285423220 | 563 41          | 285423310 | -                | -         | IVa              | 5,3                   | 50                      | 5.000 to 30.000                                      |




#### TC viscometers

- the technical measurement characteristics are in accordance with DIN 53 000, part 1, ISO 3105
- for use in combination with an automatic viscosity measuring instrument
- filling quantity: 15 to 20 ml
- overall length: ca. 355 mm
- suitable bracket Type No. 05393, Order No. 285405035

| calibrated,<br>with constant for automatic measurements |           |               |           |                |           | Capillary<br>No. | Capillary<br>Ø i [mm] | Constant<br>K<br>(approx.) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|---|-----------|---------------|-----------|----------------|-----------|------------------|-----------------------|----------------------------|--|
| Type<br>No.   | Order No. | Type<br>No.   | Order No. | Type<br>No.    | Order No. |                  |                       |                            |  |
| +10 to +80 °C   |           | -40 to +30 °C |           | +70 to +150 °C |           |                  |                       |                            |  |
| 567 03  | 285423420 | -             | -         | -              | -         | 0c               | 0,47                  | 0.003                      | 0,5 to 3   |
| 567 10  | 285423430 | 568 10        | 285423540 | 569 10         | 285423630 | I                | 0,63                  | 0,01                       | 1,2 to 10  |
| 567 13  | 285423440 | 568 13        | 285423550 | 569 13         | 285423640 | Ic               | 0,84                  | 0,03                       | 3 to 30  |
| 567 20  | 285423450 | 568 20        | 285423560 | 569 20         | 285423650 | II               | 1,13                  | 0,1                        | 10 to 100  |
| 567 23  | 285423470 | 568 23        | 285423570 | 569 23         | 285423660 | IIc              | 1,51                  | 0,3                        | 30 to 300  |
| 567 21  | 285423460 | -             | -         | -              | -         | IIa              | 1,69                  | 0,5                        | 50 to 500  |
| 567 30  | 285423480 | 568 30        | 285423580 | 569 30         | 285423670 | III              | 2,05                  | 1                          | 100 to 1.000   |
| 567 33  | 285423500 | 568 33        | 285423590 | 569 33         | 285423680 | IIIc             | 2,7                   | 3                          | 300 to 3.000   |
| 567 31  | 285423490 | -             | -         | -              | -         | IIIa             | 3,0                   | 5                          | 500 to 5.000   |
| 567 40  | 285423510 | 568 40        | 285423600 | 569 40         | 285423690 | IV               | 3,7                   | 10                         | 1.000 to 10.000                                      |
| 567 43  | 285423530 | 568 43        | 285423620 | 569 43         | 285423700 | IVc              | 4,9                   | 30                         | 3.000 to 20.000                                      |
| 567 41  | 285423520 | 568 41        | 285423610 | -              | -         | IVa              | 5,3                   | 50                         | 5.000 to 30.000                                      |

## 2.6 Micro-Ubbelohde viscometers with TC sensors



Viscometers with suspended level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behaviour. The measuring levels are marked by TC sensors. The meniscus passage is detected due to the different conductivity of the liquid phase and the gas phase. A measurement stand of the type series AVS®/S is not required. TC viscometers can be used to determine the kinematic viscosity of all liquids with Newtonian flow behaviour.

They are especially suitable for liquids that cannot be detected with other systems: opaque and/or black and/or electrically conductive measuring samples.

Due to the electrical properties of TC sensors, it is important to ensure that a suitable type is selected for the required application temperature.

### Micro TC viscometers

- the technical measurement characteristics are in accordance with DIN 53 000, Part 1
- for use in combination with an automatic viscosity measuring instrument
- filling quantity: 3 to 4 ml
- overall length: approx. 350 mm
- suitable bracket Type No. 05393, Order No. 285405035


| calibrated,<br>with constant and calibration certificate for automatic<br>measurements |           |                |           |                 |           | Capillary<br>No. | Capillary<br>Ø i [mm] | Constant K<br>(approx.) | Measuring<br>range [mm <sup>2</sup> /s]<br>(approx.) |
|--|-----------|----------------|-----------|-----------------|-----------|------------------|-----------------------|-------------------------|--|
| Type<br>No.  | Order No. | Type<br>No.    | Order No. | Type<br>No.     | Order No. |                  |                       |                         |  |
| +10 ... +80 °C   |           | -40 ... +30 °C |           | +70 ... +150 °C |           |                  |                       |                         |  |
| 572 10   | 285423710 | 573 10         | 285423780 | 574 10          | 285423850 | M I              | 0,40                  | 0,01                    | 0,4 to 6   |
| 572 13   | 285423720 | 573 13         | 285423790 | 574 13          | 285423860 | M Ic             | 0,52                  | 0,03                    | 1,2 to 18  |
| 572 20   | 285423730 | 573 20         | 285423800 | 574 20          | 285423870 | M II             | 0,70                  | 0,1                     | 4 to 60  |
| 572 23   | 285423740 | 573 23         | 285423810 | 574 23          | 285423880 | M IIc            | 0,95                  | 0,3                     | 12 to 180  |
| 572 30   | 285423750 | 573 30         | 285423820 | 574 30          | 285423890 | M III            | 1,26                  | 1                       | 40 to 800  |

## 2.7 Micro-Ubbelohde viscometers (DIN)

Viscometers with suspended level for determination of absolute and relative kinematic viscosity of liquids with Newtonian flow behavior. Due to their design, these viscometers are especially suitable for measurement of small liquid quantities

and for particularly short running times.

All viscometers are provided with ring marks. This ensures that viscometers for automatic measurements can also be checked by means of manual measurements.



### Micro-Ubbelohde viscometers (DIN)

The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18.

- in accordance with DIN 53 000, Part 1
- filling quantity: 3 to 4 ml
- overall length: approx. 290 mm

| calibrated, with constant and calibration certificate for manual measurement |           | calibrated, with constant and calibration certificate for automatic measurement |           | not calibrated, without calibration certificate; for determination of relative viscosity |           | Capillary No. | Capillary Ø i [mm] | Constant K (approx.) | Measuring range [mm <sup>2</sup> /s] (approx.) |
|--|-----------|---|-----------|--|-----------|---------------|--------------------|----------------------|--|
| Type No.   | Order No. | Type No.  | Order No. | Type No.   | Order No. |               |                    |                      |  |
| 536 10   | 285401009 | 537 10  | 285401103 | 538 10   | 285401206 | M I           | 0,40               | 0,01                 | 0,4 to 6                                       |
| 536 13   | 285401017 | 537 13  | 285401111 | 538 13   | 285401214 | M Ic          | 0,52               | 0,03                 | 1,2 to 18                                      |
| 536 11   | 285401050 | 537 11  | 285401150 | 538 11   | 285401220 | M Ia          | 0,60               | 0,05                 | 2 to 30  |
| 536 20   | 285401025 | 537 20  | 285401128 | 538 20   | 285401222 | M II          | 0,70               | 0,1                  | 4 to 60  |
| 536 23   | 285401033 | 537 23  | 285401136 | 538 23   | 285401239 | M IIc         | 0,95               | 0,3                  | 12 to 180                                      |
| 536 21   | 285401030 | 537 21  | 285401130 | 538 21   | 285401230 | M IIa         | 1,07               | 0,5                  | 20 to 300                                      |
| 536 30   | 285401041 | 537 30  | 285401144 | 538 30   | 285401247 | M III         | 1,26               | 1                    | 40 to 800                                      |

## Viscometers for dilution viscometry

Viscometers with suspended level designed according to the principle of the Ubbelohde viscometers for determination of the limit viscosity number of polymers. The limit viscosity number can be determined automatically in

combination with one of our piston burettes TITRONIC® 300 or TITRONIC® 500.

- filling quantity: 15 to 75 ml
- overall length: approx. 290 mm



| not calibrated without calibration certificate<br>Model with glass filter and discharge tube on request |           | Capillary No. | Capillary Ø i [mm] | Constant K (approx.) | Measuring range [mm <sup>2</sup> /s] (approx.) |
|---|-----------|---------------|--------------------|----------------------|--|
| Type No.  | Order No. |               |                    |                      |  |
| 531 00  | 285401403 | 0             | 0,36               | 0,001                | 0,35 to 0,6                                    |
| 531 03  | 285401428 | 0c            | 0,47               | 0,003                | 0,5 to 2                                       |
| 531 01  | 285401411 | 0a            | 0,53               | 0,005                | 0,8 to 3                                       |
| 531 10  | 285401436 | I             | 0,64               | 0,01                 | 1,2 to 6                                       |
| 531 13  | 285401444 | Ic            | 0,84               | 0,03                 | 3 to 20  |
| 531 20  | 285401452 | II            | 1,15               | 0,1                  | 10 to 60                                       |

## 2.8 Cannon-Fenske viscometers



### Cannon-Fenske routine viscometers

- comply with standards ISO 3105, BS 188 and ASTM D446 with respect to technical measuring specifications.
- are suitable for all Newtonian liquids with a viscosity of 0.35 to 20,000 mm<sup>2</sup>/s
- the present design has a deepening in the lower bend. Accordingly, these viscometers mit Abfallsystem can also be used for automatic measurements.
- filling quantity: approx. 7 to 10 ml
- overall length: approx. 245 mm

| calibrated,<br>with constant and<br>calibration certificate<br>for manual measurements |           | calibrated,<br>with constant and calibration<br>certificate<br>for automatic measurements |           | Capillary<br>No. | Capillary<br>Ø i [mm] | Constant K<br>(Richtwert) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|--|-----------|---|-----------|------------------|-----------------------|---------------------------|--|
| Type No.   | Order No. | Type No.  | Order No. |                  |                       |                           |  |
| 513 00   | 285403507 | 520 00  | 285403704 | 25               | 0,30                  | 0,002                     | 0,4 to 1,6   |
| 513 03   | 285403515 | 520 03  | 285403712 | 50               | 0,44                  | 0,004                     | 0,8 to 3,2   |
| 513 01   | 285403523 | 520 01  | 285403729 | 75               | 0,54                  | 0,008                     | 1,6 to 6,4   |
| 513 10   | 285403531 | 520 10  | 285403737 | 100              | 0,63                  | 0,015                     | 3 to 15  |
| 513 13   | 285403548 | 520 13  | 285403745 | 150              | 0,78                  | 0,035                     | 7 to 35  |
| 513 20   | 285403556 | 520 20  | 285403753 | 200              | 1,01                  | 0,1                       | 20 to 100  |
| 513 23   | 285403564 | 520 23  | 285403761 | 300              | 1,27                  | 0,25                      | 50 to 200  |
| 513 21   | 285403572 | 520 21  | 285403778 | 350              | 1,52                  | 0,5                       | 100 to 500   |
| 513 30   | 285403589 | 520 30  | 285403786 | 400              | 1,92                  | 1,2                       | 240 to 1.200   |
| 513 33   | 285403597 | 520 33  | 285403794 | 450              | 2,35                  | 2,5                       | 500 to 2.500   |
| 513 40   | 285403601 | 520 40  | 285403807 | 500              | 3,20                  | 8                         | 1.600 to 8.000                                       |
| 513 43   | 285403618 | 520 43  | 285403815 | 600              | 4,20                  | 20                        | 4.000 to 20.000                                      |



### Cannon-Fenske reverse flow viscometers


- Comply with standards ISO 3105, ASTM D 446 with respect to technical measuring specifications.
- filling quantity: approx. 12 ml
- overall length: approx. 295 mm

| calibrated,<br>with 3 ring marks,<br>with 2 constants and calibration certificate,<br>only for manual measurement |           | Capillary No. | Capillary<br>Ø i [mm] | Constant K<br>(approx.) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|---|-----------|---------------|-----------------------|-------------------------|--|
| Type No.  | Order No. |               |                       |                         |  |
| 511 00  | 285403001 | 25            | 0,31                  | 0,002                   | 0,4 to 1,6   |
| 511 03  | 285403018 | 50            | 0,42                  | 0,004                   | 0,8 to 3,2   |
| 511 01  | 285403026 | 75            | 0,54                  | 0,008                   | 1,6 to 6,4   |
| 511 10  | 285403034 | 100           | 0,63                  | 0,015                   | 3 to 15  |
| 511 13  | 285403042 | 150           | 0,78                  | 0,035                   | 7 to 35  |
| 511 20  | 285403059 | 200           | 1,02                  | 0,1                     | 20 to 100  |
| 511 23  | 285403067 | 300           | 1,26                  | 0,25                    | 50 to 200  |
| 511 21  | 285403075 | 350           | 1,48                  | 0,5                     | 100 to 500   |
| 511 30  | 285403083 | 400           | 1,88                  | 1,2                     | 240 to 1200  |
| 511 33  | 285403091 | 450           | 2,20                  | 2,5                     | 500 to 2500  |
| 511 40  | 285403104 | 500           | 3,10                  | 8                       | 1.600 to 8.000                                       |
| 511 43  | 285403112 | 600           | 4,00                  | 20                      | 4.000 to 20.000                                      |

Cannon-Fenske routine viscometers with additional tube and threads

comply with standards ISO 3105, BS 188 with respect to technical measuring specifications. These viscometers are preferably used for automatic measurements. The additional filling and cleaning tube and the glass thread ensure safe operational use. The calibrated viscometers are delivered with manufacturer's certificate in accordance with DIN 55 350, Part 18.

- are suitable for all Newtonian liquids with a viscosity of 0.35 to 20,000 mm<sup>2</sup>/s.
- filling quantity: approx. 7 to 12 ml
- overall length: approx. 245 mm



| calibrated,<br>with constant and calibration certificate for<br>automatic measurements |           | Capillary No. | Capillary<br>Ø i [mm] | Constant K<br>(approx.) | Measuring range<br>[mm <sup>2</sup> /s]<br>(approx.) |
|--|-----------|---------------|-----------------------|-------------------------|--|
| Type No.   | Order No. |               |                       |                         |  |
| 546 00   | 285402116 | 25            | 0,30                  | 0,002                   | 0,4 to 1,6   |
| 546 03   | 285402132 | 50            | 0,44                  | 0,004                   | 0,8 to 3,2   |
| 546 01   | 285402124 | 75            | 0,54                  | 0,008                   | 1,6 to 6,4   |
| 546 10   | 285402149 | 100           | 0,63                  | 0,015                   | 3 to 15  |
| 546 13   | 285402157 | 150           | 0,78                  | 0,035                   | 7 to 35  |
| 546 20   | 285402165 | 200           | 1,01                  | 0,1                     | 20 to 100  |
| 546 23   | 285402181 | 300           | 1,27                  | 0,25                    | 50 to 200  |
| 546 21   | 285402173 | 350           | 1,52                  | 0,5                     | 100 to 500   |
| 546 30   | 285402198 | 400           | 1,92                  | 1,2                     | 240 to 1.200   |
| 546 33   | 285402202 | 450           | 2,35                  | 2,5                     | 500 to 2.500   |
| 546 40   | 285402219 | 500           | 3,20                  | 8                       | 1.600 to 8.000                                       |
| 546 43   | 285402227 | 600           | 4,20                  | 20                      | 4.000 to 20.000                                      |

## 2.9 Ostwald viscometers



### Ostwald viscometers

- filling quantity: 3 ml
- overall length: approx. 220 mm
- only available without calibration

| with ring marks,<br>without calibration certificate,<br>for manual measurements |           | Capillary<br>Ø i [mm] | Transit time for water<br>~ [s] | Constant K<br>(approx.) | for use from<br>[mm <sup>2</sup> /s]<br>(approx.) |
|---|-----------|-----------------------|---------------------------------|-------------------------|---|
| Type No.  | Order No. |                       |                                 |                         |   |
| 509 03  | 285404006 | 0,3                   | 250                             | 0,004                   | 0,3   |
| 509 04  | 285404014 | 0,4                   | 75                              | 0,01                    | 1   |
| 509 05  | 285404022 | 0,5                   | 30                              | 0,03                    | 2,5   |
| 509 06  | 285404039 | 0,6                   | 15                              | 0,07                    | 5,5   |
| 509 07  | 285404047 | 0,7                   | 10                              | 0,1                     | 10  |



### Micro-Ostwald viscometers

- are suitable for measurements of small liquid quantities even with tendency to excessive foam formation.
- filling quantity: 2 ml
- overall length: approx. 290 mm

| calibrated,<br>with ring marks,<br>with constant and<br>calibration certificate<br>for manual<br>measurements |           | calibrated,<br>with ring marks,<br>with constant and<br>calibration certificate<br>for automatic<br>measurements |           | not calibrated,<br>without calibration<br>certificate;<br>for determination of<br>relative viscosity |           | Capillary<br>No. | Capillary<br>Ø i [mm] | Constant K<br>(approx.) | Measuring<br>range [mm <sup>2</sup> /s]<br>(approx.) |
|---|-----------|--|-----------|--|-----------|------------------|-----------------------|-------------------------|--|
| Type No.  | Order No. | Type No.   | Order No. | Type No.   | Order No. |                  |                       |                         |  |
| 516 10  | 285404203 | 517 10   | 285404306 | 518 10   | 285404409 | I                | 0,43                  | 0,01                    | 0,4 bis 6  |
| 516 13  | 285404211 | 517 13   | 285404314 | 518 13   | 285404417 | Ic               | 0,60                  | 0,03                    | 1,2 bis 18   |
| 516 20  | 285404228 | 517 20   | 285404322 | 518 20   | 285404425 | II               | 0,77                  | 0,1                     | 4 bis 60   |
| 516 23  | 285404236 | 517 23   | 285404339 | 518 23   | 285404433 | IIc              | 1,00                  | 0,3                     | 12 bis 180   |
| 516 30  | 285404244 | 517 30   | 285404347 | 518 30   | 285404441 | III              | 1,36                  | 1                       | 40 bis 800   |

## 3 Accessories

### 3.1.1 Viscometers - Brackets and stands

#### Brackets and stands

All brackets and stands are designed to ensure that the viscometers are held vertically. They also protect the viscometers from breakage. The maximum deviation is  $< 1^\circ$ .

In conjunction with Xylem Analytics Germany and other commercially available seethrough thermostats the viscometers can only be used with the appropriate stand or bracket.

For DIN Ubbelohde viscometers that are used as reference measuring standard, specifically modified bracket (VZ 5840, see below) must be used.

#### Brackets made of stainless steel suitable for use with all Ubbelohde viscometers

for manual and automatic measurements

| Type No.   | Order No. |
|--|-----------|
| 053 92   | 285405043 |
| VZ 5840 (accessory for reference measuring standard) | 285417201 |

#### suitable for use with Ubbelohde viscometers with TC sensors

| Type No. | Order No. |
|----------|-----------|
| 053 93   | 285405035 |

#### suitable for use with Cannon-Fenske reverse flow viscometers

for manual measurements (not illustrated)

| Type No. | Order No. |
|----------|-----------|
| 053 96   | 285405019 |

#### suitable for use with Micro-Ostwald viscometers

for manual and automatic measurements

| Type No. | Order No. |
|----------|-----------|
| 053 97   | 285405027 |



## Accessories for brackets for reference measuring standard

DIN Ubbelohde viscometers which are used as testing standard


should be stored in a specially modified viscometer bracket according to official inspection / calibration authorities. The extension set for the test standard (VZ 5840) as supplement to the bracket

guarantees vertical slope with a maximum deviation of  $< 1^\circ$  and the centered positioning of the capillaries.

| Type No. | Order No. |
|----------|-----------|
| VZ 5840  | 285417201 |



Centering plate and clamp for viscometers



Holder made of polyamid



### Polyamide bracket

for use with Cannon-Fenske routine viscometers, Cannon-Fenske reverse flow viscometers and all Ostwald viscometers for manual measurements only

| Typ-Nr. | Bestell-Nr. |
|---------|-------------|
| 064 99  | 285405105   |

### PTFE bracket

for use with Cannon-Fenske routine viscometers, for automatic measurements only (not illustrated), to be applied in measuring stands AVS®/S, AVS®/SK


| Typ-Nr. | Bestell-Nr. |
|---------|-------------|
| 065 99  | 285405113   |

## 3.1.2 Temperature stabilization jackets

### Temperature stabilization jackets

In the absence of a see-through thermostat the temperature of capillary viscometers can be stabilized in this type of jacket using circulation thermostats in the temperature range 0 to 180 °C.

The support plate has been designed to facilitate changing the viscometer when required. An additional hole is provided in the support plate so that a control thermometer can be fitted. A quick-action seal simplifies changing viscometers.



### Temperature stabilization jacket with support plate for Ubbelohde viscometers

| Type No.               | Order No. | Item No. | Comment   |
|------------------------|-----------|----------|---|
| 577 00                 | 285405508 |          | complete, without viscometer  |
| <b>Component parts</b> |           |          |   |
| 577 01                 | 285405516 | 1        | temperature stabilization jacket, straight  |
| 238 00                 | 285424130 | 2        | support plate with 4 silicone rings for Ubbelohde viscometers (d = 4, 6, 8 and 10 mm) |
| 225 34                 | 285405532 | 3        | silicone O-ring, ND 60  |
| 072 34                 | 285405549 | 4        | quick-action seal, NW 60  |

### Silicone rings

| Type No. | Order No. | d mm | D mm | h mm |
|----------|-----------|------|------|------|
| 228 11   | 285405808 | 4    | 10   | 5    |
| 228 14   | 285405816 | 6    | 16   | 5    |
| 228 16   | 285405824 | 8    | 16   | 5    |
| 228 17   | 285405832 | 10   | 16   | 5    |

## 3.1.3 LabPump

### LabPump

The LabPump VZ 5655 (not illustrated) used with manual and semi-automatic measurements to extract and pump solutions:

- For manual measurements: To suck sample to the measuring bulb.
- For semi-automatic AVS®-Systems: To discharge sample to a waste bottle.

Since the LabPump VZ 5655 and the connections are made of PTFE or stainless steel, the pump is suitable for use with aggressive mediums.

The discharging of the sample from the viscometer is possible for viscosities up to 30000 mm<sup>2</sup>/s. Additional to the LabPump, the discharge set type no VZ 5624 is required.

| Type No.      | Order No. |
|---------------|-----------|
| VZ 5655, 230V | 1040755   |
| VZ 5665, 115V | 1040757   |
| VZ 5624       | 285414845 |

## 3.2 A. for Measurement Devices

### 3.2.1 Control thermometers

#### Control thermometers

Precision thermometers type TFX 430 of brand Ebro® with high accuracy  $\pm 0.05$  °C to control bath temperature. The thermometer is equipped with an adapter for installation to bath thermostats of series CT72 and CT52.

| Type No. | Order No. | Article             | Description  |
|----------|-----------|---------------------|--|
| VZ 7330  | 285421110 | Control thermometer | Pt100 resistance thermometer, with manufacturer's calibration  |
| VZ 7340  | 285421120 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 20 °C  |
| VZ 7341  | 285421130 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 25 °C  |
| VZ 7342  | 285421140 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 30 °C  |
| VZ 7343  | 285421170 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 40 °C  |
| VZ 7344  | 285421180 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 100 °C |
| VZ 7345  | 285421290 | Control thermometer | Pt100 resistance thermometer, with DAkkS calibration at 135 °C |

## Technical Data

- Measuring range: -100 to +500 °C
- Accuracy:  $\pm 0.05$  °C from -50 to +199.99 °C
- Operating temperature: -20 °C to +50 °C
- Sensor length: 200 mm
- Four-wire system
- Water proof acc. IP 67
- Battery operated, battery lifetime approx. 5 years
- Battery exchangeable

The available control thermometers only differ regarding the calibration: Manufacturer's calibration (VZ 7330) respective DAkkS calibration for certain temperatures, which are commonly used in viscometry.



## 3.2.2 Safety devices

### Safety bottle for suction mode

For operation in suction mode, AVS® systems have to be protected by a safety bottle. The bottle, being positioned between capillary tube of viscometer and ViscoPump, avoids the entering of liquid into the control unit: In case of malfunction, the sample is first hauled into the viscometer.

The safety bottle without glass screws, VZ 7022, for use with silicon caps and hoses, is contained in the supplement set VZ 8526. By using the enclosed hose, the standard hose set VZ 5505 for pressure mode can be switched to suction mode.

Beside of the safety bottle VZ 7022, there is also safety bottle VZ 7021 with glass screws, to connect with flanged PTFE hoses and screw cap connectors. The PTFE hoses are mainly used for aggressive samples as sulfuric acid, which would attack silicon.

With PTF tubes there is - depending on the application - a broad variety due to different tube length, diameter and screw connectors, so that they cannot be described in this catalogue. Please contact for special requests your sales partner for SI Analytics® viscometry systems.



| Type No. | Order No. | Description   |
|----------|-----------|---|
| VZ 7022  | 285420277 | Safety bottle without glass screws                                    |
| VZ 8526  | 285420530 | Upgrading set "suction" for ViscoPump in addition to hose set VZ 5505 |
| VZ 7021  | 285420269 | Safety bottle with glass screws                                       |

## Safety sensors

When using the safety bottle VZ 7021 / VZ 7022, the safety device for vacuum line VZ 8552 is strongly recommended. This capacitive sensor is mounted under the safety bottle and - in case of malfunction - registers liquid, resulting in a shutdown of the ViscoPump.

In case of an attached waste system, additionally the safety sensor VZ 8551 can be used: This safety sensor detects the weight of the waste bottle and, in case of overfill, switches off the pumping of the waste liquid.



| Type No. | Order No. | Description                                       |
|----------|-----------|---|
| VZ 8552  | 1054303   | Safety sensor for suction line, capacitive sensor |
| VZ 8551  | 1054112   | Safety sensor for waste bottle                    |



## 3.2.3 Filtration system ProClean II



### Filtration system ProClean II

For filtration of solutions in viscometry. All standard solutions can be filtrated, also with aggressive solvents used in polymer analytics. The system is configured for 100 ml screw neck bottles GL 45. The solutions are filtrated by using wire-mesh filters made of Hastelloy or filter disks made of vitreous PE, fixed on a holder by a screw cap. The sample solutions are filtrated by sucking through the wire mesh by using disposable syringes with luer-fitting (e.g. 20ml, not in scope of supply). As filtration is done by suction, there is no risk of splashing - important for filtration of dangerous liquids. After use, the filter blanks can be cleaned and reused.

The single components of the filtration system can be ordered separately.



Reusable filter VZ 7094, VZ 7095



Disposable filter VZ 7097, VZ 7098

| Type No | Order No. | Article                              | Description  |
|---------|-----------|--------------------------------------|--|
| VZ 7090 | 285422500 | Filtration system ProClean II, Set-5 | 5x sample bottles with screw caps and filter holders; filter blanks (Hastelloy) 20µm and 30µm mesh size (10pcs. per set) |
| VZ 7092 | 285422470 | Filter holder for ProClean II        | 5 pcs.   |
| VZ 7093 | 285422510 | Screw caps for ProClean II           | For fixing of filters, 10 pcs  |
| VZ 7094 | 285422480 | Filter 20µm mesh size                | Hastelloy C4, 10 pcs   |
| VZ 7095 | 285422520 | Filter 30µm mesh size                | Hastelloy C4, 10 pcs   |
| VZ 7096 | 285422490 | Bottle set for ProClean II           | 5x sample bottles with screw caps (with and without bore hole)   |
| VZ 7097 | 285422590 | Filter disks for ProClean II         | PE, 50 pcs   |
| VZ 7098 | 285422600 | Filter disks for ProClean II         | PE, 500W pcs   |

## 3.2.4 Thermostat vessel for ViscoClock plus

The standards for capillary viscometry require an accuracy for temperature control of  $\pm 0.02$  K (ISO 3104, DIN 53 000, ASTM D445) respective  $\pm 0.05$  for relative measurements (ISO 1628-1)


To achieve this accurate tempering, for viscometry we recommend our high-precision glass-panelled bath thermostats.

As an alternative, in case of the ViscoClock plus we offer a special double wall vessel for tempering. This vessel made of DURAN® glass is connected to a thermostat bath circulator.

The temperature control vessel can also be used for manual viscosity measurement, using a stopwatch and our viscometer brackets (e.g. 053 92).

Due to unavoidable heat loss and for safety reasons, the range of measuring temperature is restricted to  $\pm 15$  °C compared to room temperature.

For better temperature control, the use of an additional magnetic stirrer is recommended.



### Specifications

- Temperature range: +10 to + 40 °C
- Temperature accuracy:  $\pm 0.05$  K
- Diameter of tubing olive: 10.5 mm
- Recommended: Stirring using magnetic stirrer, length of stirring bar approx. 40 mm.

| Type No. | Order No. | Description                                    |
|----------|-----------|--|
| VZ 6574  | 285424400 | Temperature control vessel for ViscoClock plus |

## 3.2.5 Burette TITRONIC® 500 for sample preparation

For viscometry of polymer solutions, sample preparation has to be done first. As an alternative to volumetric flasks, SI Analytics® offers the piston burette TITRONIC® 500 together with an exchangeable head. For dosing aggressive or high-viscous solvents, e.g. sulfuric acid or phenol/dichlorobenzene, a special and robust exchangeable head WA 50 V was designed.

There are the following advantages for sample preparation with the TITRONIC® 500 in comparison to using volumetric flasks:

- No manual work with aggressive solvents.
- When preparing in a volumetric flask, the weighing has to be done precisely to reach a target value, e.g. 250 mg when using a 50 ml flask for a concentration of 0.5 g/dl. The precise weighing is difficult for polymer samples, e.g. granules. With TITRONIC® 500, the dosed solvent volume matches exactly to the sample weight - therefore it becomes obsolete to reach exact target sample weight.
- When using TITRONIC® 500, standard sample bottles can be used instead of volumetric flasks - no more filling-up to the graduation mark.
- In case that a stirring bar is necessary for dissolution, in volumetric flasks it has to be taken out and adhering liquid has to be rinsed to the flask. This is not necessary with the piston burette, as target volume is dosed before dissolving.
- If sample is heated for dissolution, in volumetric flasks the solution needs to cool down to room temperature before fill-up.
- High reproducibility ( $\pm 0.2\%$ ) of dosing volume.
- Simple and inexpensive compared to systems using gravimetric solvent dosing.
- Interchangeable unit WA 50 V with highly resistant dosing piston, ensuring tightness also for aggressive and high viscous solvents (e.g. 96 % sulfuric acid).




Beside of the sample weight, in the software of TITRONIC® 500 also the content of foreign substances, e.g. glass fibers, can be typed in. This percentage of sample weight is disregarded for calculation of the amount of solvent. Therefore the solvent amount is calculated on the pure polymer weight, according to viscometry standards, e.g. ISO 1628 or ISO 307.



| Type No. | Order No. | Article  | Description  |
|----------|-----------|--|--|
| T 500    | 285220200 | TITRONIC® 500  | TITRONIC® 500 basic unit without magnetic stirrer TM 235, with stand rod and titration clamp Z 305, and power supply 100-240 V |
| T 500-M1 | 285220210 | TITRONIC® 500 with magnetic stirrer  | TITRONIC® 500 basic unit with magnetic stirrer TM 235, with stand rod and titration clamp Z 305, and power supply 100-240 V    |
| WA 50 V  | 285220360 | 50 ml interchangeable unit for aggressive solvents, e.g. conc. sulfuric acid                       | with brown glass bottle for titrant, GL 45 and S 40-bottle adapter, tubes, drip tube and titration tip                         |
| WA 50    | 285220350 | 50 ml interchangeable unit for aqueous solutions or non-aggressive solvents, e.g. ethanol, acetone | with brown glass bottle for titrant, GL 45 and S 40-bottle adapter, tubes, drip tube and titration tip                         |

## 3.2.6 AVS® measuring stands



### AVS® measuring stands

Measuring stands of the series AVS®/S can be used to measure the flow-through time in viscometers automatically. The measuring stands can be connected to all measuring instruments made by Xylem Analytics Germany for automatic measurement of viscosity and operate with all standard viscometers for repetitive measurements (except Ostwald viscometer type series 509 xx).

Automatic measurements have the following advantages:

- the repetitive standard deviation is less than for manual measurements
- the measurement is free from subjective factors of influence
- the results can be printed and/or be automatically documented with a data memory system
- automatic processing of sample series is available.

The measuring stands or brackets can be exchanged, having only negligible influence to the flow time:

The distance between the levels of the automatic optoelectronic unloading system is 40.00 mm  $\pm$  0.03 mm. Due to the low tolerance, the exchange of a measuring stand results in a deviation of only  $\pm$  0.05 % for Ubbelohde viscometers.

For repetitive measurements with viscosity measuring instruments and Ubbelohde viscometers with measuring stands, the standard deviation ca.  $\pm$  0.03 % for stable and particle-free sample liquids.

Manually calibrated Ubbelohde viscometers can also be used in AVS® measuring stands. If the automatic sensing levels do not correspond to the ring marks, the height difference of the meniscus detection system will result in a changed viscometer constant. The difference amounts to 0.1 % per millimeter of height offset for Ubbelohde viscometers.

| Type No.                   | Order No. | Description  |
|----------------------------|-----------|--|
| Measuring stand AVS®/S     | 285410502 | Metal measuring stand AVS®/S, preferably for nonaqueous bath fluids                            |
| Measuring stand AVS®/SK    | 285410876 | PVDF measuring stand AVS®/SK, corrosion-free, suitable for aqueous and nonaqueous bath fluids  |
| Measuring stand AVS®/SK-CF | 285410892 | PVDF measuring stand AVS®/SK-CF, particularly for the use of Cannon-Fenske routine viscometers |
| Measuring stand AVS®/SK-V  | 285410905 | PVDF measuring stand AVS®/SK-V, particularly for the use of dilution viscometers               |

## Measuring stands

|   | AVS®/S   | AVS®/S-HT      | AVS®/SK               | AVS®/SK-CF                            | AVS®/SK-V  |
|---|--|----------------|-----------------------|---------------------------------------|--|
| Available viscometers                     | Ubbelohde viscometers in accordance with DIN, ASTM, ISO 3105, Micro-Ubbelohde viscometers, Micro-Ostwald viscometers |                |                       | Cannon-Fenske-routine viscometer      | Ubbelohde-dilution viscometer  |
| Temperature range                         | -80 to +80 °C  | -80 to +200 °C | -40 to +80 °C         | -80 to +80 °C                         | -40 to +80 °C other temperature ranges available on request                              |
| Suitable brackets (type no.)              | 05392<br>05397   |                |                       | no bracket required                   |  |
| Material                                  | Aluminium, TiO <sub>2</sub> -anodized  |                | PVDF, stainless steel | Aluminium, TiO <sub>2</sub> -anodized | PVDF, stainless steel  |
| Dimensions (W x H x D) mm                 | 90 x 447 x 90  | 90 x 496 x 90  | 90 x 447 x 90         | 90 x 447 x 90                         | 90 x 447 x 90  |
| Weight (kg) appr.                         | 1.0  | 1.25           | 0.8                   | 1.0                                   | 0.8  |
| Accessories included in scope of delivery | Bracket Type No. 05392 for Ubbelohde viscometers, tube/cable combination VZ 5505                                     |                |                       | tube/cable combination VZ 5505        | tube/cable combination VZ 5857, magnetic stirring rods, fastening springs for viscometer |

**Note:**

When TC viscometers are being used, a bracket type no. 053 93, with the necessary tube set is required only. A measuring stand is not required.


Suitable for use with the measuring units: AVS® 370, AVS® 470, AVS®Pro (optical detection)

Suitable for use with the thermostatic baths: CT 72/P, CT 72/2-TT, CT 72/2, CT 72/4

Electrical connection: Cable VZ 6225 for all measuring stands to all instruments (is included in hose sets VZ 5505, VZ 5622 and VZ 5857), control lamp as function display

Distance between measuring levels: 40.00 mm ± 0.03 mm at 25 °C

Signal transmission: Optically using optical fibres from the measuring level in the stand head, converted into analogue signal from stand to measuring instrument



# More SI Analytics Products

## Titration

### Dosing, titrating and water analysis according to Karl Fischer can be so easy

#### The Titrators of the TitroLine® series

Based on our know-how, we have developed a range of new reliable laboratory instruments for dosing, titrating and Karl Fischer water analysis: The burettes TITRONIC® 300 and 500 and the titrators TitroLine® 5000, 7000, 7500 KF, 7500 KF trace and the universal titrators TitroLine® 7750 and 7800. These instruments combine easy handling with maximum accuracy, and the robustness required for the daily operation in the laboratory.

- ▶ For the completely successful laboratory operation we also offer a wide range of accessories perfectly supporting the titrators with all their functions, such as the sample changers TW alpha plus and TW 7400, the burette TITRONIC® 500.



## ▲ The new Titrator TitroLine® 7800 gives even more options

The TitroLine® 7800 is as well as the TitroLine® 7750 the all-rounder for both potentiometric titration and volumetric KF titration. Compared with the TitroLine® 7750 the TitroLine® 7800 offers the opportunity to also connect digital IDS sensors.



All-rounder: TitroLine® 7800



OptiLine 6

### OptiLine 6 sensor

Many titration applications and methods, e.g. N Ph.Eur or USP prescribe the use of an indicator for the titration end point. There are also methods that explicitly require the use of a photometric sensor. The OptiLine 6 is a new photometric sensor that can be used like any other sensor. Thanks to the additional analog BNC / DIN connection, it can be connected to any titrator or even a pH meter with an appropriate measuring input.

# pH measurement

## Perfectly matching measurement systems

### Benchtop and handheld meters

We offer several product lines for a wide range of applications. From the beginner class to the digital multi-channel device you get everything from a single source.



*HandyLab MKII*



*HandyLab 845  
Ideal for the entry to pH measurement*



*Multifunctional ProLab 2500 and 5000*



### Laboratory electrodes

Suitable for our versatile devices, we offer countless types of electrodes for almost every application. These electrodes are still produced by hand at our location in Mainz in order to ensure the best possible quality and a long shelf life and precise measurement.



## Process electrodes, armatures and accessories

The reliable measuring of pH, ORP, conductivity and D.O. values up to temperature within the process requires individual solutions. Our extensive range of process electrodes includes all applications for measurements in aqueous solutions in the temperature range from  $-30\text{ }^{\circ}\text{C}$  to  $140\text{ }^{\circ}\text{C}$  at a pressure up to 12 bar. Furthermore, many of our electrodes are registered according to the ATEX guidelines 94/9/EG. The retractable holders and their control system enable flexible measurements with an ideal positioning of the electrode in the medium.



*SteamLine Electrode*



*CHEMtrac 810*

## Hotplates and Stirrers

The laboratory hotplates from SI Analytics have the benefits of the glass-ceramic heating surface which has proven by millions in households. Chemical resistance, a high-grade surface quality and a resistance to temperature shock of more than  $700\text{ }^{\circ}\text{C}$  provide the user maximum benefits compared to conventional hot plate materials. The always plane and pore-free surface enables even most stubborn dirt to be removed.



*Stirrer SLR*





# Xylem |'zīləm|

- 1) The tissue in plants that brings water upward from the roots;
- 2) a leading global water solutions company.


Xylem is a leading global water solutions company dedicated to advancing sustainable impact and empowering the people who make water work every day. Xylem connects diverse capabilities and innovative technologies to provide tailored solutions across the entire water cycle. From moving, treating and measuring water to optimizing and maintaining water systems, Xylem collaborates with customers to solve their most critical challenges. Together, through partnerships with utilities, industrial manufacturers, building operators, and communities, we are building a more water-secure world.

For more information on how Xylem can help you, go to [www.xylem.com](http://www.xylem.com)



Xylem Analytics Germany Sales GmbH & Co. KG  
Erich-Dombrowski-Str. 4  
55127 Mainz  
Germany

Tel +49 6131 894-5111  
Info.XAGS@xylem.com  
[www.xylemanalytics.com](http://www.xylemanalytics.com)



Please find the contact details of our sales, service and other departments on our website:  
[www.xylemanalytics.com/en/contact/general-contact](http://www.xylemanalytics.com/en/contact/general-contact)

SI Analytics® is a registered trademark of Xylem Inc. or one of its subsidiaries.  
Technical changes reserved. Prices are subject to change  
© 2026 Xylem Analytics Germany Sales GmbH & Co. KG.  
980 087US

Version 10/2019

